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- Dr. Sara Stead, Central Science Laboratory (CSL), Sand Hutton, York YO41 1LZ, England ;
- Dr. Steve Lehotay, USDA Agricultural Research Service (ARS), Eastern Regional Research Center, Microbial Biophysics and Residue Chemistry Research Unit, 600 East Mermaid Lane, Wyndmoor, PA 19038, USA.

Rapporteurs :

- Dr. Joe Boison, Canadian Food Inspection Agency, 116 Veterinary Road, S7N 2R3 Canada ;
- Dr. Michael Murphy, University of Minnesota.

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11. REFERENCES BIBLIOGRAPHIQUES

European Commission. 2002. Commission Decision (EC) N° 2002/657 of 12 August 2002 implementing Council Directive 96/23/EC concerning the performance of analytical methods and interpretation of results. Official Journal of European Communities. L221: 8-36.

Bahmani K, Shahbazi Y, Nikousefat Z. 2020. Monitoring and risk assessment of tetracycline residues in foods of animal origin. Food Science and Biotechnology. 29:441-448.

Berendsen BJA, Pikkemaat MG, Stolker LAM. 2011. Are antibiotic screening approaches sufficiently adequate? A proficiency test. Anal. Chim. Acta. 685:170-175.

Cantwell H, O'Keefe M. 2006. Evaluation of the Premi Test and comparison with the One-Plate Test for the detection of antimicrobials in kidney. Food Additives and Contaminants. 23:120-125.

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Commission E. 2002. Commission Decision (EC) N° 2002/657 of 12 August 2002 implementing Council Directive 96/23/EC concerning the performance of analytical methods and interpretation of results. Off J Eur Comm. L221: 8-36.

CRL. 2010. Guideline for the validation of screening methods for residues of veterinary medicines (initial validation and transfer). Available from: <http://ec.europa.eu/food/food/chemicalsafety/residues/lab_analysis_en.htm>: Guideline_Validation_Screening_en.pdf. 1-18.

Dang PK, Degand G, Danyi S, Pierret G, Delahaut P, Ton VD, Maghuin-Rogister G, Scippo M-L. 2010. Validation of a two-plate microbiological method for screening antibiotic residues in shrimp tissue. Anal. Chim. Acta. 672:30-39.

Dang PK, Degand G, Douny C, Ton VD, Maghuin-Rogister G, Scippo ML. 2011. Optimisation of a new two-plate screening method for the detection of antibiotic residues in meat. Int. J. Food Sci. Technol. 46:2070-2076.

Ekli R, Adzitey F, Huda N. 2019. Prevalence of resistant Salmonella spp. isolated from raw meat and liver of cattle in the Wa Municipality of Ghana. IOP Conference Series: Earth and Environmental Science. 287:012006.

Ezenduka EV, Ike OS, Anaelom NJ. 2014. Rapid detection of antimicrobial residues in poultry: A consequence of non-prudent use of antimicrobials. Health. 6:149-152.

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Ezenduka EV, Okorie-Kanu OJ, Nwanta JA. 2019. Comparative analysis of two microbiological tests in the detection of oxytetracycline residue in chicken using ELISA as gold standard. *Journal of Immunoassay and Immunochemistry*. 40:617-629.

Gaudin V, Gaugain-Juhel M, Moretain JP, Sanders P. 2008. AFNOR validation of Premi®Test, a microbiological-based screening tube-test for the detection of antimicrobial residues in animal muscle tissue. *Food Addit. & Contam.: Part A*. 25:1451 - 1464.

Gaudin V, Hedou C, Rault A, Sanders P, Verdon E. 2009. Comparative study of three screening tests, two microbiological tube tests, and a multi-sulphonamide ELISA kit for the detection of antimicrobial and sulphonamide residues in eggs. *Food Addit. & Contam.: Part A*. 26:427–440.

Gondová Z, Kožárová I, Poláková Z, Mad'arová M. 2014. Comparison of Four Microbiological Inhibition Tests for the Screening of Antimicrobial Residues in the Tissues of Food-Producing Animals. *Italian Journal of Animal Science*. 13:3521.

Hoff R, Ribarcki F, Zancanaro I, Castellano L, Spier C, Barreto F, Fonseca SH. 2012. Bioactivity-based screening methods for antibiotics residues: A comparative study of commercial and inhouse developed kits. *Food Additives and Contaminants - Part A Chemistry, Analysis, Control, Exposure and Risk Assessment*. 29:577-586.

Kilinc B, Cakli S. 2008. Screening for antibiotic residues in the trout by the Four Plate test, Premi test and ELISA test. *Eur. Food Res. Technol*. 226:795-799.

Kilinc B, Meyer C, Hilge V. 2007. Evaluation of the EEC four-plate test and Premi test for screening antibiotic residues in trout (*Salmo trutta*). *Int. J. Food Sci. Technol*. 42:625-628.

Kim DP, Degand G, Douny C, Pierret G, Delahaut P, Ton VD, Granier B, Scippo M-L. 2013. Preliminary Evaluation of Antimicrobial Residue Levels in Marketed Pork and Chicken Meat in the Red River Delta Region of Vietnam. *Food and Public Health*. 3 267-276.

Kirrella GAK, Deeb AMM, Abdallah RMI. 2017. Safety of frozen liver for human consumption. *Journal of Food and Drug Analysis*. 25:520-524.

Kozarova I, Janosova J, Mate D, Tkacikova S. 2009. Evaluation of three different microbial inhibition tests for the detection of sulphamethazine residues in the edible tissues of rabbit. *Food Additives and Contaminants*. Vol. 26,:978-987.

Kožárová I, Juščáková D, Šimková J, Milkovičová M, Kožár M. 2020. Effective screening of antibiotic and coccidiostat residues in food of animal origin by reliable broad-spectrum residue screening tests. *Italian Journal of Animal Science*. 19:487-501.

Kožárová I, Marcinčák S, Reitznerová A, Bartkovský M, Mačanga J, Marcinčáková D, Klempová T, Čertík M, *Verification for the presence of inhibitory substances in poultry meat after the consumption of the feed mixture supplemented with fermented feed*. 2016.

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Mitiku B, Adujna M. 2018. Antimicrobial residue occurrence and its public health risk of beef meat in Debre Tabor and Bahir Dar, Northwest Ethiopia. *Veterinary World.* 11:902-908.

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Nhung NT, Van NTB, Cuong NV, Duong TTQ, Nhat TT, Hang TTT, Nhi NTH, Kiet BT, Hien VB, Ngoc PT, Campbell J, Thwaites G, Carrique-Mas J. 2018. Antimicrobial residues and resistance against critically important antimicrobials in non-typhoidal *Salmonella* from meat sold at wet markets and supermarkets in Vietnam. *Int. J. Food Microbiol.* 266:301-309.

Njoga EO, Onunkwo JI, Okoli CE, Ugwuoke WI, Nwanta JA, Chah KF. 2018. Assessment of antimicrobial drug administration and antimicrobial residues in food animals in Enugu State, Nigeria. *Tropical Animal Health and Production.* 50:897-902.

Omeiza GK, Ajayi IE, Ode OJ, *Assessment of antimicrobial drug residues in beef in Abuja, the Federal Capital Territory, Nigeria.* 2012; Vol. 48, p 283-9.

Pikkemaat MG, Rapallini MLBA, Oostra-van Dijk S, Elferink JWA. 2009. Comparison of three microbial screening methods for antibiotics using routine monitoring samples. *Anal. Chim. Acta.* 637:298-304.

Pikkemaat MG, Rapallini MLBA, Zuidema T, Elferink JWA, Oostra-van Dijk S, Driessen-van Lankveld WDM. 2011. Screening methods for the detection of antibiotic residues in slaughter animals: comparison of the European Union Four-Plate Test, the Nouws Antibiotic Test and the PremiTest (applied to muscle and kidney). *Food Additives and Contaminants.* 28:26-34.

Rahimi Z, Shahbazi Y, Ahmadi F. 2018. Comparative Screening of Chloramphenicol Residue in Chicken Tissues Using Four Plate Test and Premi®Test Methods. *Pharm Sci.* 24:157-162.

Rahman MR, Lou Z, Antora R, Saqib M. 2015. Lateral Flow Test Strip Approaches For Rapid Detection Of Food Contaminants And Pathogens: Review. *American Research Thoughts.* 1:1222-1248.

Rakotoharinome M, Pognon D, Randriamparany T, Ming JC, Idoumbin J-P, Cardinale E, Porphyre V. 2014. Prevalence of antimicrobial residues in pork meat in Madagascar. *Tropical Animal Health and Production.* 46:49-55.

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Sallam KI, Saad FSS, Abdelkhalek A. 2022. Health risk assessment of antimicrobial residues in sheep carcasses marketed in Kuwait. Food Chem. 383:132401.

Schneider M, Lehotay S. 2008. A comparison of the FAST, Premi® and KIS tests for screening antibiotic residues in beef kidney juice and serum. Anal. Bioanal. Chem. 390:1775-1779.

Schneider MJ, Mastovska K, Lehotay SJ, Lightfield AR, Kinsella B, Shultz CE. 2009. Comparison of screening methods for antibiotics in beef kidney juice and serum. Anal. Chim. Acta. 637:290–297.

Staruch L, Mati M. 2012. Detection of inhibitory substances in poultry meat and poultry meat quality AND POULTRY MEAT QUALIT. FOLIA VETERINARIA, 56, Supplementum I: 51—53.

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Wu Q, Shabbir MAB, Peng D, Yuan Z, Wang Y. 2021. Microbiological inhibition-based method for screening and identifying of antibiotic residues in milk, chicken egg and honey. Food Chem. 363:130074.

Yang S-C, Yu M-C, Lee Y-H, Wang J-L, *Antibiotic Residues in Meat and Eggs in Taiwan: A Local Surveillance*. 2016; Vol. 12, p 1-6.

Ziani K, Pérez-López M, Mansouri A, Khaled MB, Rodriguez AS, Slimani M. 2018. Assessment of Oxytetracycline Residue in Cooked and Raw Meat of Chicken Broilers Before and After the End of Official Withdrawal Period. Food Analytical Methods. 11:2528-2537.

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
Premi@Test

Muscle Screening test for the presence of antibiotic and sulfamide residues in particular beef meat. NF VALIDATION certification only possible with pasteurized meat (in the exception of ground meat). The following antibiotic residues were investigated: Penicillins (PCP = 4 µg/g); Amoxicillins (ACP = 10 µg/g); Chloramphenicol (CCP = 100 µg/g); Tetracyclines (TCM = 10 µg/g); Sulphonamides (SCP = 10 µg/g); Streptomycin (STP = 100 µg/g); Tylosin (TYP = 100 µg/g); Macrolides (MCP = 100 µg/g); Trimethoprim (TTP = 10 µg/g). The CCP values are determined by gelatin diffusion method. The detection limit is 0.1 mg/kg. Colistin, gentamicin, ampicillin, amphotericin, quinolones, tetracyclines, and lincosamides are detected at the maximum residual level (MRL). <http://www.biopharm.com>

Contiene: 25 ampolline con Bacillus anthracis/anthracis (B. anthracis) in a solid agar medium, dispensable gelatin plate, perforated coverlid + syringe.

Perfect description: Premi@Test is a solid agar medium intended for microbiological detection of antibiotic residues in animal products. It is a solid agar medium, dispensable gelatin plate, perforated coverlid + syringe.

Test principle: Premi@Test is based on the inhibition of the growth of Bacillus anthracis (B. anthracis) by the action of antibiotic residues. A detection limit of 0.1 mg/kg is achieved.

Caution: This test is suitable to detect antibiotic and other antibiotics, any combination with sulfamide should be avoided at all times. It is advised to wear masks throughout the whole procedure. Use clean paper or clean towel to dry hands.

- Instructions for use:**
- Wash hands thoroughly before starting the test procedure.
 - Read the requirements of application before starting the application of the adjustment solution.
 - Remove the cap carefully from the required amount of ampoules. Do not touch the ampoules or needles.
 - Take approximately 1 cm³ of the adjustment solution and add it to the ampoules. Do not touch the ampoules or needles.
 - Allow the ampoules to sit for the syringe for each sample.
 - Pipette out 10 µl of the sample into the syringe. Do not touch the ampoules or needles.
 - Insert the ampoules into the syringe and push the plunger down to the 10 µl mark. Do not touch the ampoules or needles.
 - Close the test ampoules with the cap and push the plunger down to the 10 µl mark.
 - Incubate the test ampoules in a Premi@Test incubator or water bath at 37 °C ± 0.5 °C.
 - It is recommended to use, in parallel, a negative control of the same matrix that is being tested. Read the test results after the negative control has changed color.
 - All necessary equipment (test paper, syringe, incubator, and solvent) are available in the Premi@Test Starter Kit.
 - To find out the use of the Premi@Test there is a thorough document available on request for all users.

Reading test results:

- Read the color only from the bottom 2 µl of the ampoule.
- A colorless change (from purple towards white) indicates the absence of antibiotic / sulfamide above the detection limit.
- No color change indicates the presence of antibiotic and / or sulfamide at or above the level of detection of the test.

Application: The use of an ampoule control is required (obligatory for NF VALIDATION). Start to check the color of the sample control after an incubation time of 2 h at 37 °C and then in intervals of 10 min until the color of the sample control stopped changing from purple to white. If that time, the sample should be compared to the negative control. Because a change from purple to white can lead to incorrect results of the detection capability (CCP). If there is no color change of the sample control after 4 h, repeat the test. The negative control can be prepared by heating and storing muscle from a sample of the negative control. New use water is a negative control.

Positive control: It is strongly recommended to test a positive control (control) is available from Biopharm AG to verify the correct application of the test.

Storage: The ampoules should be stored correctly (30 °C ± 2 °C, Attention: DO NOT FREEZE).


Label liability: Premi@Test is a screening test and such 100% accuracy of the test results cannot be guaranteed. Besides, the assessment of color, in particular that of a yellow (purple) color, may differ from person to person, in cases where severe color changes are involved for the user, test results should be confirmed by a suitable complementary analysis method. Biopharm makes no warranty of any kind, either expressed or implied, except that the materials from which the products are made are of standard quality.

Recommendations for sample preparation: of fish, shrimp, egg, honey, etc. urine, blood and feces are available at Biopharm AG. These matrices require special adaptation by the customer. Support is available at Biopharm upon request.

Biopharm AG
 Postfach 101111, Postfach 101111
 An der Ammersee 30, Postfach 101111
 84529 Gars am Inn, Germany
 Tel.: +49 (0) 84 41 - 93 10 0
 Fax: +49 (0) 84 41 - 93 10 10
 E-Mail: info@biopharm.com
www.biopharm.com

Autobrief/Supermarket Board:
 Dr. Ralf W. Ziemer (Vorsitzender/Chairman)
 Vorstand/Supermarket Board
 Chairman, Dr. Hans-Peter, Julia Henrich, Dr. Peter Schabert
 Anhangsleiter/Committee Members:
 Anhangsleiter/Committee Members
 UStD Dr. W. W. Dr. W. Dr. W. Dr. W.

Premi@Test

Muscle Screening test for the presence of antibiotic and sulfamide residues in particular beef meat. NF VALIDATION certification only possible with pasteurized meat (in the exception of ground meat). The following antibiotic residues were investigated: Penicillins (PCP = 4 µg/g); Amoxicillins (ACP = 10 µg/g); Chloramphenicol (CCP = 100 µg/g); Tetracyclines (TCM = 10 µg/g); Sulphonamides (SCP = 10 µg/g); Streptomycin (STP = 100 µg/g); Tylosin (TYP = 100 µg/g); Macrolides (MCP = 100 µg/g); Trimethoprim (TTP = 10 µg/g). The CCP values are determined by gelatin diffusion method. The detection limit is 0.1 mg/kg. Colistin, gentamicin, ampicillin, amphotericin, quinolones, tetracyclines, and lincosamides are detected at the maximum residual level (MRL). <http://www.biopharm.com>

Contiene: 25 ampolline con Bacillus anthracis/anthracis (B. anthracis) in a solid agar medium, dispensable gelatin plate, perforated coverlid + syringe.

Perfect description: Premi@Test is a solid agar medium intended for microbiological detection of antibiotic residues in animal products. It is a solid agar medium, dispensable gelatin plate, perforated coverlid + syringe.

Test principle: Premi@Test is based on the inhibition of the growth of Bacillus anthracis (B. anthracis) by the action of antibiotic residues. A detection limit of 0.1 mg/kg is achieved.

Caution: This test is suitable to detect antibiotic and other antibiotics, any combination with sulfamide should be avoided at all times. It is advised to wear masks throughout the whole procedure. Use clean paper or clean towel to dry hands.

- Instructions for use:**
- Wash hands thoroughly before starting the test procedure.
 - Read the requirements of application before starting the application of the adjustment solution.
 - Remove the cap carefully from the required amount of ampoules. Do not touch the ampoules or needles.
 - Take approximately 1 cm³ of the adjustment solution and add it to the ampoules. Do not touch the ampoules or needles.
 - Allow the ampoules to sit for the syringe for each sample.
 - Pipette out 10 µl of the sample into the syringe. Do not touch the ampoules or needles.
 - Insert the ampoules into the syringe and push the plunger down to the 10 µl mark. Do not touch the ampoules or needles.
 - Close the test ampoules with the cap and push the plunger down to the 10 µl mark.
 - Incubate the test ampoules in a Premi@Test incubator or water bath at 37 °C ± 0.5 °C.
 - It is recommended to use, in parallel, a negative control of the same matrix that is being tested. Read the test results after the negative control has changed color.
 - All necessary equipment (test paper, syringe, incubator, and solvent) are available in the Premi@Test Starter Kit.
 - To find out the use of the Premi@Test there is a thorough document available on request for all users.

Reading test results:

- Read the color only from the bottom 2 µl of the ampoule.
- A colorless change (from purple towards white) indicates the absence of antibiotic / sulfamide above the detection limit.
- No color change indicates the presence of antibiotic and / or sulfamide at or above the level of detection of the test.

Application: The use of an ampoule control is required (obligatory for NF VALIDATION). Start to check the color of the sample control after an incubation time of 2 h at 37 °C and then in intervals of 10 min until the color of the sample control stopped changing from purple to white. If that time, the sample should be compared to the negative control. Because a change from purple to white can lead to incorrect results of the detection capability (CCP). If there is no color change of the sample control after 4 h, repeat the test. The negative control can be prepared by heating and storing muscle from a sample of the negative control. New use water is a negative control.

Positive control: It is strongly recommended to test a positive control (control) is available from Biopharm AG to verify the correct application of the test.

Storage: The ampoules should be stored correctly (30 °C ± 2 °C, Attention: DO NOT FREEZE).


Label liability: Premi@Test is a screening test and such 100% accuracy of the test results cannot be guaranteed. Besides, the assessment of color, in particular that of a yellow (purple) color, may differ from person to person, in cases where severe color changes are involved for the user, test results should be confirmed by a suitable complementary analysis method. Biopharm makes no warranty of any kind, either expressed or implied, except that the materials from which the products are made are of standard quality.

Recommendations for sample preparation: of fish, shrimp, egg, honey, etc. urine, blood and feces are available at Biopharm AG. These matrices require special adaptation by the customer. Support is available at Biopharm upon request.

Biopharm AG
 Postfach 101111, Postfach 101111
 An der Ammersee 30, Postfach 101111
 84529 Gars am Inn, Germany
 Tel.: +49 (0) 84 41 - 93 10 0
 Fax: +49 (0) 84 41 - 93 10 10
 E-Mail: info@biopharm.com
www.biopharm.com

Autobrief/Supermarket Board:
 Dr. Ralf W. Ziemer (Vorsitzender/Chairman)
 Vorstand/Supermarket Board
 Chairman, Dr. Hans-Peter, Julia Henrich, Dr. Peter Schabert
 Anhangsleiter/Committee Members:
 Anhangsleiter/Committee Members
 UStD Dr. W. W. Dr. W. Dr. W. Dr. W.

Premi@Test

Muscle Screening test for the presence of antibiotic and sulfamide residues in particular beef meat. NF VALIDATION certification only possible with pasteurized meat (in the exception of ground meat). The following antibiotic residues were investigated: Penicillins (PCP = 4 µg/g); Amoxicillins (ACP = 10 µg/g); Chloramphenicol (CCP = 100 µg/g); Tetracyclines (TCM = 10 µg/g); Sulphonamides (SCP = 10 µg/g); Streptomycin (STP = 100 µg/g); Tylosin (TYP = 100 µg/g); Macrolides (MCP = 100 µg/g); Trimethoprim (TTP = 10 µg/g). The CCP values are determined by gelatin diffusion method. The detection limit is 0.1 mg/kg. Colistin, gentamicin, ampicillin, amphotericin, quinolones, tetracyclines, and lincosamides are detected at the maximum residual level (MRL). <http://www.biopharm.com>

Contiene: 25 ampolline con Bacillus anthracis/anthracis (B. anthracis) in a solid agar medium, dispensable gelatin plate, perforated coverlid + syringe.

Perfect description: Premi@Test is a solid agar medium intended for microbiological detection of antibiotic residues in animal products. It is a solid agar medium, dispensable gelatin plate, perforated coverlid + syringe.

Test principle: Premi@Test is based on the inhibition of the growth of Bacillus anthracis (B. anthracis) by the action of antibiotic residues. A detection limit of 0.1 mg/kg is achieved.

Caution: This test is suitable to detect antibiotic and other antibiotics, any combination with sulfamide should be avoided at all times. It is advised to wear masks throughout the whole procedure. Use clean paper or clean towel to dry hands.

- Instructions for use:**
- Wash hands thoroughly before starting the test procedure.
 - Read the requirements of application before starting the application of the adjustment solution.
 - Remove the cap carefully from the required amount of ampoules. Do not touch the ampoules or needles.
 - Take approximately 1 cm³ of the adjustment solution and add it to the ampoules. Do not touch the ampoules or needles.
 - Allow the ampoules to sit for the syringe for each sample.
 - Pipette out 10 µl of the sample into the syringe. Do not touch the ampoules or needles.
 - Insert the ampoules into the syringe and push the plunger down to the 10 µl mark. Do not touch the ampoules or needles.
 - Close the test ampoules with the cap and push the plunger down to the 10 µl mark.
 - Incubate the test ampoules in a Premi@Test incubator or water bath at 37 °C ± 0.5 °C.
 - It is recommended to use, in parallel, a negative control of the same matrix that is being tested. Read the test results after the negative control has changed color.
 - All necessary equipment (test paper, syringe, incubator, and solvent) are available in the Premi@Test Starter Kit.
 - To find out the use of the Premi@Test there is a thorough document available on request for all users.

Reading test results:

- Read the color only from the bottom 2 µl of the ampoule.
- A colorless change (from purple towards white) indicates the absence of antibiotic / sulfamide above the detection limit.
- No color change indicates the presence of antibiotic and / or sulfamide at or above the level of detection of the test.

Application: The use of an ampoule control is required (obligatory for NF VALIDATION). Start to check the color of the sample control after an incubation time of 2 h at 37 °C and then in intervals of 10 min until the color of the sample control stopped changing from purple to white. If that time, the sample should be compared to the negative control. Because a change from purple to white can lead to incorrect results of the detection capability (CCP). If there is no color change of the sample control after 4 h, repeat the test. The negative control can be prepared by heating and storing muscle from a sample of the negative control. New use water is a negative control.

Positive control: It is strongly recommended to test a positive control (control) is available from Biopharm AG to verify the correct application of the test.

Storage: The ampoules should be stored correctly (30 °C ± 2 °C, Attention: DO NOT FREEZE).


Label liability: Premi@Test is a screening test and such 100% accuracy of the test results cannot be guaranteed. Besides, the assessment of color, in particular that of a yellow (purple) color, may differ from person to person, in cases where severe color changes are involved for the user, test results should be confirmed by a suitable complementary analysis method. Biopharm makes no warranty of any kind, either expressed or implied, except that the materials from which the products are made are of standard quality.

Recommendations for sample preparation: of fish, shrimp, egg, honey, etc. urine, blood and feces are available at Biopharm AG. These matrices require special adaptation by the customer. Support is available at Biopharm upon request.

Biopharm AG
 Postfach 101111, Postfach 101111
 An der Ammersee 30, Postfach 101111
 84529 Gars am Inn, Germany
 Tel.: +49 (0) 84 41 - 93 10 0
 Fax: +49 (0) 84 41 - 93 10 10
 E-Mail: info@biopharm.com
www.biopharm.com

Autobrief/Supermarket Board:
 Dr. Ralf W. Ziemer (Vorsitzender/Chairman)
 Vorstand/Supermarket Board
 Chairman, Dr. Hans-Peter, Julia Henrich, Dr. Peter Schabert
 Anhangsleiter/Committee Members:
 Anhangsleiter/Committee Members
 UStD Dr. W. W. Dr. W. Dr. W. Dr. W.

Premi@Test

Muscle Screening test for the presence of antibiotic and sulfamide residues in particular beef meat. NF VALIDATION certification only possible with pasteurized meat (in the exception of ground meat). The following antibiotic residues were investigated: Penicillins (PCP = 4 µg/g); Amoxicillins (ACP = 10 µg/g); Chloramphenicol (CCP = 100 µg/g); Tetracyclines (TCM = 10 µg/g); Sulphonamides (SCP = 10 µg/g); Streptomycin (STP = 100 µg/g); Tylosin (TYP = 100 µg/g); Macrolides (MCP = 100 µg/g); Trimethoprim (TTP = 10 µg/g). The CCP values are determined by gelatin diffusion method. The detection limit is 0.1 mg/kg. Colistin, gentamicin, ampicillin, amphotericin, quinolones, tetracyclines, and lincosamides are detected at the maximum residual level (MRL). <http://www.biopharm.com>

Contiene: 25 ampolline con Bacillus anthracis/anthracis (B. anthracis) in a solid agar medium, dispensable gelatin plate, perforated coverlid + syringe.

Perfect description: Premi@Test is a solid agar medium intended for microbiological detection of antibiotic residues in animal products. It is a solid agar medium, dispensable gelatin plate, perforated coverlid + syringe.

Test principle: Premi@Test is based on the inhibition of the growth of Bacillus anthracis (B. anthracis) by the action of antibiotic residues. A detection limit of 0.1 mg/kg is achieved.

Caution: This test is suitable to detect antibiotic and other antibiotics, any combination with sulfamide should be avoided at all times. It is advised to wear masks throughout the whole procedure. Use clean paper or clean towel to dry hands.

- Instructions for use:**
- Wash hands thoroughly before starting the test procedure.
 - Read the requirements of application before starting the application of the adjustment solution.
 - Remove the cap carefully from the required amount of ampoules. Do not touch the ampoules or needles.
 - Take approximately 1 cm³ of the adjustment solution and add it to the ampoules. Do not touch the ampoules or needles.
 - Allow the ampoules to sit for the syringe for each sample.
 - Pipette out 10 µl of the sample into the syringe. Do not touch the ampoules or needles.
 - Insert the ampoules into the syringe and push the plunger down to the 10 µl mark. Do not touch the ampoules or needles.
 - Close the test ampoules with the cap and push the plunger down to the 10 µl mark.
 - Incubate the test ampoules in a Premi@Test incubator or water bath at 37 °C ± 0.5 °C.
 - It is recommended to use, in parallel, a negative control of the same matrix that is being tested. Read the test results after the negative control has changed color.
 - All necessary equipment (test paper, syringe, incubator, and solvent) are available in the Premi@Test Starter Kit.
 - To find out the use of the Premi@Test there is a thorough document available on request for all users.

Reading test results:

- Read the color only from the bottom 2 µl of the ampoule.
- A colorless change (from purple towards white) indicates the absence of antibiotic / sulfamide above the detection limit.
- No color change indicates the presence of antibiotic and / or sulfamide at or above the level of detection of the test.

Application: The use of an ampoule control is required (obligatory for NF VALIDATION). Start to check the color of the sample control after an incubation time of 2 h at 37 °C and then in intervals of 10 min until the color of the sample control stopped changing from purple to white. If that time, the sample should be compared to the negative control. Because a change from purple to white can lead to incorrect results of the detection capability (CCP). If there is no color change of the sample control after 4 h, repeat the test. The negative control can be prepared by heating and storing muscle from a sample of the negative control. New use water is a negative control.

Positive control: It is strongly recommended to test a positive control (control) is available from Biopharm AG to verify the correct application of the test.

Storage: The ampoules should be stored correctly (30 °C ± 2 °C, Attention: DO NOT FREEZE).

Label liability: Premi@Test is a screening test and such 100% accuracy of the test results cannot be guaranteed. Besides, the assessment of color, in particular that of a yellow (purple) color, may differ from person to person, in cases where severe color changes are involved for the user, test results should be confirmed by a suitable complementary analysis method. Biopharm makes no warranty of any kind, either expressed or implied, except that the materials from which the products are made are of standard quality.

Recommendations for sample preparation: of fish, shrimp, egg, honey, etc. urine, blood and feces are available at Biopharm AG. These matrices require special adaptation by the customer. Support is available at Biopharm upon request.

Biopharm AG
 Postfach 101111, Postfach 101111
 An der Ammersee 30, Postfach 101111
 84529 Gars am Inn, Germany
 Tel.: +49 (0) 84 41 - 93 10 0
 Fax: +49 (0) 84 41 - 93 10 10
 E-Mail: info@biopharm.com
www.biopharm.com

Autobrief/Supermarket Board:
 Dr. Ralf W. Ziemer (Vorsitzender/Chairman)
 Vorstand/Supermarket Board
 Chairman, Dr. Hans-Peter, Julia Henrich, Dr. Peter Schabert
 Anhangsleiter/Committee Members:
 Anhangsleiter/Committee Members
 UStD Dr. W. W. Dr. W. Dr. W. Dr. W.

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
Fast and reliable test for the detection of antibiotic residues and sulfonamides particularly in dairy cows. The certification NF VALIDATION proves that the results are reliable, accurate and validated in dairy cows. The antibiotics concerned are the following:

- Penicillins: Benzathine (CIP = 500 µg/kg), ampicillin (CIP = 100 µg/kg), ampicillin (CIP = 100 µg/kg)
- Tetracyclines: Chloramphenicol (CIP = 50 µg/kg), doxycycline (CIP = 50 µg/kg), oxytetracycline (CIP = 100 µg/kg)
- Sulfonamides: Sulfadiazine (CIP = 100 µg/kg), sulfamonomethyl (CIP = 100 µg/kg)
- Metronidazole: Metronidazole (CIP = 100 µg/kg)

The value of CIP is the abbreviation of the content of the units of sale. The detection capacities are other residues of veterinary medicines. Chloramphenicol, tetracyclines, quinolones, furazolidone, florfenicol and trimethoprim are classified in the same way as the others.

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Premi®Test

Microbial Screening test for the presumptive detection of antibiotic residues in meat, fish, shrimps, eggs, kidneys, liver, urine, blood and feed. NF VALIDATION certification scope includes beef, pork and poultry meat (with the exception of ground meat).

RBP 31/02 - 04/11 ALTERNATIVE ANALYTICAL METHODS FOR AGRIBUSINESS
Certified by AFNOR Certification
End of validity: 30th August 2018
<http://nf-validation.afnor.org>

2014/10/15 - insert 1011

Contents

25 ampoules with *Bacillus stearothermophilus* var. *calidolactis* in a solid agar medium, disposable pipette tips, perforated cover foil + syringe.

Product description

Premi®Test is a broad spectrum microbial screening test especially developed for the detection of antimicrobial substances, such as antibiotic and sulphonamide residues, in fresh meat, fish and eggs at or below most Maximum Residue Level (MRL).

Test principle

Premi®Test is based on the inhibition of the growth of *Bacillus stearothermophilus*, a micro-organism very sensitive to many antibiotics and sulphonamide residues.

A standardised number of spores are imbedded in an agar medium with selected nutrients. When meat fluid is added to the Premi®Test and heated at 64 °C, the spores will germinate. The germinated spores will multiply and form an acid when no inhibitory substances are present. This will be visible by a colour change from purple to yellow of the indicator in the ampoule. When antimicrobial residues are present in sufficient amount (above the detection level) no growth will occur and the colour will remain purple.

Special "sample-procedures" for fish, shrimps, eggs, kidneys, liver, urine, blood and feed are available at R-Biopharm AG.

Caution

This test is extremely sensitive to antibiotics and other inhibitory substances, any contamination with such materials should be avoided at all time. It is advised to wash hands thoroughly before starting the test procedure. Use tissue paper or a clean towel to dry hands.

Instruction for use

- Wash hands thoroughly before starting the test procedure.
- Cut out the required number of ampoules without damaging the aluminum foil of the adjacent ampoules.
- Remove foil carefully from the required amount of ampoules (do not open more ampoules as needed).
- Take approximately 2 cm³ of lean meat and use a meat press to extract approximately 250µl of meat fluid. Alternatively, the meat juice can be obtained with the help of the Multipress or by freeze/thawing the meat. (A technical bulletin on the Multipress is available at www.r-biopharm.com).
- Use a new disposable tip on the syringe for each sample.
- Pipette 100 µl of the fluid onto the agar in the ampoule. Do not distort the agar.
- Allow to stand at room temperature for 20 minutes for a pre-incubation.
- Flush the meat juice away by gently filling and emptying the test ampoule twice with demineralized (demi) water. **Wash the ampoules with demineralized water only and do not use tap-water!**
- Remove the last water carefully from the test ampoule.
- Close the test ampoule with the foil supplied to avoid evaporation.
- Incubate the test ampoule in a Premi®Test Incubator or water bath (64 °C ± 0,5 °C).
- It is **strongly recommended** to use, in parallel, a negative control of the same matrix that is being tested. Read the test results after the negative control has changed color.
- All necessary equipment (meat press, scissors, incubator, and timer) are available in the Premi®Test Starter Kit.
- To facilitate the use of the Premi®Test there is a Powerpoint document available on request for all users.

Reading test results

- Read the color only from the bottom 2/3 part of the ampoule.
- A clear color change (from purple towards yellow) indicates the absence of antibiotics / sulphonamides above the detection limit.
- No clear color change indicates the presence of antibiotics and/or sulphonamides at/or above the limit of detection of the test.

Negative control

The use of a negative control is **strongly recommended (obligatory for NF VALIDATION)**. Start to check the color of the negative control after an incubation time of 2 h 40 min and then in intervals of 5 min until the color of the negative control has changed from purple to yellow. If there is no color change of the negative control after 3 h 10 min, repeat the test. The negative control can be prepared by freezing and storing meat juice from a sample tested negative. **Never use water as a negative control!**

Positive control

It is recommended to test regularly a positive control (protocol is available from R-Biopharm AG) to verify the correct application of the test.

Storage

The ampoules should be stored cool (3 – 10 °C). Attention: **DO NOT FREEZE!**

Limited Liability

Premi®Test is a screening test and as such 100% accuracy of the test results cannot be guaranteed. Besides, the assessment of colour, in particular that of a yellow/purple result, may differ from person to person. In cases where severe consequences are involved for the user, test results should be confirmed by a validated comprehensive analytical method. The data corresponds to our present state of technology and provides information on our products and their uses. R-Biopharm makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. Defective products will be replaced. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. R-Biopharm shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

R-Biopharm AG
Postanschrift / Postal Address:
An der neuen Bergstraße 17
64297 Darmstadt, Germany
Sitz / Corporate Seat: Pfungstadt
Tel.: +49 (0) 61 51 - 81 02-0
Fax: +49 (0) 61 51 - 81 02-40
E-mail: Info@r-biopharm.de
www.r-biopharm.com

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Premi®Test

RBP 31/02-04/11
ALTERNATIVE ANALYTICAL
METHODS FOR AGRIBUSINESS
Certified by AFNOR Certification
End of validity: 30th August 2022
<http://nf-validation.afnor.org>

Microbial Screening test for the presumptive detection of antibiotic and sulfonamide residues in particularly fresh meat. NF VALIDATION certification scope includes beef, pork and poultry meat (with the exception of ground meat). The following selected antibiotics were investigated:
Penicillins: Penicillin G (CC β = 6 μ g/kg); Amoxicillin (CC β = 11 μ g/kg); Cloxacillin (CC β = 150 μ g/kg)
Tetracyclines: Chlorotetracycline (CC β = 160 μ g/kg); Oxytetracycline (CC β = 160 μ g/kg); Doxycycline (CC β = 100 μ g/kg)
Sulfonamides: Sulfadimethoxine (CC β = 75 μ g/kg); Sulfadiazine (CC β = 90 μ g/kg)
Macrolides: Erythromycine A (CC β = 200 μ g/kg); Tylosine A (CC β = 90 μ g/kg)
The CC β values were determined by spiking in meat juice; the detection capabilities in other samples may differ.
Cephalosporins, Aminoglycosides, Amphenicols, Quinolones, Tiamulin, and Lincomycin are detected above the Maximum Residue Level (MRL).
2019/03/21 - Insert 1011

Contents

25 ampoules with *Bacillus stearothermophilus* var. *calidolactis* in a solid agar medium, disposable pipette tips, perforated cover foil + syringe.

Product description

Premi®Test is a broad spectrum microbial screening test especially developed for the detection of antimicrobial substances, such as antibiotic and sulfonamide residues particularly in fresh meat, in the range of MRL.

Test principle

Premi®Test is based on the inhibition of the growth of *Bacillus stearothermophilus*, a micro-organism very sensitive to many antibiotic and sulfonamide residues.

A standardised number of spores are imbedded in an agar medium with selected nutrients. When meat fluid is added to the Premi®Test and heated at 64 °C, the spores will germinate. The germinated spores will multiply and form an acid when no inhibitory substances are present. This will be visible by a colour change from purple to yellow of the indicator in the ampoule. When antimicrobial residues are present in sufficient amount (above the detection level) no growth will occur and the colour will remain purple.

Caution

This test is sensitive to certain antibiotics and other inhibitory substances, any contamination with such materials should be avoided at all time. It is advised to wash hands thoroughly before starting the test procedure. Use tissue paper or a clean towel to dry hands.

Instruction for use

- Wash hands thoroughly before starting the test procedure.
- Cut out the required number of ampoules without damaging the aluminum foil of the adjacent ampoules.
- Remove foil carefully from the required amount of ampoules (do not open more ampoules as needed).
- Take approximately 2 cm³ of lean meat and use a meat press to extract approximately 250 μ l of meat fluid. Alternatively, the meat juice can be obtained with the help of the Multipress or by freeze/thawing the meat. (A technical bulletin on the Multipress is available at www.r-biopharm.com).
- Use a new disposable tip on the syringe for each sample.
- Pipette 100 μ l of the fluid onto the agar in the ampoule. Do not distort the agar.
- Allow to stand at room temperature for 20 minutes for a pre-incubation.
- Flush the meat juice away by gently filling and emptying the test ampoule twice with demineralized (demi) water. Wash the ampoules with demineralized water only and do not use tap-water.
- Remove the last water carefully from the test ampoule.
- Close the test ampoule with the foil supplied to avoid evaporation.
- Incubate the test ampoule in a Premi®Test Incubator or water bath (64 °C \pm 1 °C).
- It is required to use, in parallel, a negative control of the same matrix that is being tested. Read the test results after the negative control has changed color.
- All necessary equipment (meat press, scissors, incubator, and timer) are available in the Premi®Test Starter Kit.
- To facilitate the use of the Premi®Test there is a Powerpoint document available on request for all users.

Reading test results

- Read the color only from the bottom 2/3 part of the ampoule.
- A clear color change (from purple towards yellow) indicates the absence of antibiotics / sulfonamides above the detection limit.
- No clear color change indicates the presence of antibiotics and/or sulfonamides at/or above the limit of detection of the test.

Negative control

The use of a negative control is required (obligatory for NF VALIDATION). Start to check the color of the negative control after an incubation time of 2 h 40 min and then in intervals of 5 min until the color of the negative control has changed from purple to yellow. At that time, the samples should be compared to the negative control, because a delayed measuring (> 5 min.) can lead to increased values of the detection capability (CC β). If there is no color change of the negative control after 4 h, repeat the test. The negative control can be prepared by freezing and storing meat juice from a sample tested negative. Never use water as a negative control!

Positive control

It is strongly recommended to test a positive control (protocol is available from R-Biopharm AG) to verify the correct application of the test.

Storage

The ampoules should be stored cool (3 – 10 °C). Attention: DO NOT FREEZE!

Limited Liability

Premi®Test is a screening test and as such 100% accuracy of the test results cannot be guaranteed. Besides, the assessment of colour, in particular that of a yellow/purple result, may differ from person to person. In cases where severe consequences are involved for the user, test results should be confirmed by a validated comprehensive analytical method. R-Biopharm makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality.

Recommendations for sample-preparations of fish, shrimps, eggs, kidney, liver, urine, blood and feed are available at R-Biopharm AG. These matrices require individual validation by the customer. Supportive information are available at R-Biopharm upon request.

R-Biopharm AG
Postanschrift / Postal Address:
An der neuen Bergstraße 17
64297 Darmstadt, Germany
Sitz / Corporate Seat: Pfungstadt
Tel.: +49 (0) 61 51 - 81 02-0
Fax: +49 (0) 61 51 - 81 02-40
E-mail: info@r-biopharm.de
www.r-biopharm.com

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Annexe 3. Certificat No. 060601 d'Octobre 2021.

CERTIFICATION
AOAC® Performance TestedSM


Certificate No.
060601

The AOAC Research Institute hereby certifies the performance of the method known as:

Premi®Test

manufactured by
DSM PREMI® Test B. V.
P. O. Box 6500
6401 JH HEERLEN
The Netherlands

This method has been evaluated in the AOAC® Performance Tested MethodsSM Program and found to perform as stated by the manufacturer contingent to the comments contained in the manuscript. This certificate means that an AOAC® Certification Mark License Agreement has been executed which authorizes the manufacturer to display the AOAC Performance TestedSM certification mark along with the statement - "THIS METHOD'S PERFORMANCE WAS REVIEWED BY AOAC RESEARCH INSTITUTE AND WAS FOUND TO PERFORM TO THE MANUFACTURER'S SPECIFICATIONS" - on the above-mentioned method for a period of one calendar year from the date of this certificate (October 23, 2021 – December 31, 2022). Renewal may be granted at the end of one year under the rules stated in the licensing agreement.

 _____ Scott Coates, Senior Director Signature for AOAC Research Institute	October 23, 2021 _____ Date
--	-----------------------------------

2275 Research Blvd., Ste. 300, Rockville, Maryland, USA Telephone: +1-301-924-7077 Fax: +1-301-924-7089
Internet e-mail: aoacri@aoac.org * World Wide Web Site: <http://www.aoac.org>

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METHOD AUTHORS Angelique de Rijk	MANUFACTURER DSM Nutritional Products DSM Premi® Test BV P. O. Box 6500 6401 JH Heerlen The Netherlands	CERTIFICATION MARK LICENSÉE R-Biopharm AG An der Neuen Bergstraße 17, 64297 Darmstadt, Germany
METHOD NAME Premi® Test	CATALOG NUMBERS R3925, R3900	
INDEPENDENT LABORATORY Central Science Laboratory (CSL) York YO41 1LZ, England USDA Agricultural Research Service Eastern Regional Research Center Microbial Biophysics and Residue Chemistry Research Unit 600 East Mermaid Lane Wyndmoor, PA 19038, USA	AOAC EXPERTS AND PEER REVIEWERS Joe Boison ¹ , Michael Murphy ² ¹ Canadian Food Inspection Agency, Saskatoon, SK, Canada ² University of Minnesota, Minneapolis, MN, USA	
APPLICABILITY OF METHOD Target analyte – Penicillin G Matrix – Bovine muscle tissue Performance claims – See Table 1		
ORIGINAL CERTIFICATION DATE September 27, 2006	CERTIFICATION RENEWAL RECORD Renewed annually through December 2022.	
METHOD MODIFICATION RECORD 1. November 2019 Level 1	SUMMARY OF MODIFICATION 1. Editorial changes to inserts.	
Under this AOAC® Performance Tested™ License Number, 060601 this method is distributed by: NONE	Under this AOAC® Performance Tested™ License Number, 060601 this method is distributed as: NONE	

PRINCIPLE OF THE METHOD (1)

The Premi®test is based on the inhibition of the growth of *Bacillus Stearothermophilus*, a bacterium which is very sensitive to many antibiotics and sulfa compounds. A standardized number of spores are imbedded in an agar medium with selected nutrients. A small amount of meat juice is added to the test tube. After a pre-incubation period of 20 minutes at room temperature (18-22 °C i.e. 64.4-71.6 °F), the test tube is washed and incubated for approximately 3 hours at 64°C (147°F). At this temperature, the spores will start to grow (germinate) and multiply with the production of acid if no antimicrobial (non-inhibitory) substances are present in the meat juice. This will be visible by a color change in the test tube from purple to yellow. When antimicrobial (inhibitory) substances are present (at or above the detection sensitivity), no growth will occur and the color in the test tube will remain purple.

DISCUSSION OF THE VALIDATION STUDY (1)

The results show that PremiTest met all the requirements for this AOAC validation study for the detection of penicillin G residues in bovine muscle tissue. The results also show that there is cross reactivity with many other antimicrobial compounds. PremiTest is a broad-spectrum screening test which can detect many relevant antimicrobials in meat.

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R-biopharm Premi®Test, AOAC® Performance Tested™ Certification Number 060601

Table 1: Summary of performance parameters of Premi®Test for penicillin G obtained in the study (1)

Laboratory	Study	Performance Pen G (ppb)	Remarks
CSL	Part 1: 6 out of 6	10	A range of 5 different concentrations from 0 – 50 ppb was tested. At 10 ppb, all 6 of 6 samples analysed tested positive.
USDA	Part 2: 30 out of 30	10	Results confirmed that the detection sensitivity of the PremiTest for pen G was 10 ppb with 95% confidence level
DSM	Part 3: Ruggedness, interference, cross-reactivity	10	The presence of other antimicrobial drugs did not interfere with the PremiTest kits ability to detect pen G residues at a concentration of 10 ppb added to bovine muscle tissue.
USDA	Part 4: Incurred tissue	10	The detection sensitivity determined using fortified tissue samples was further confirmed using incurred bovine muscle tissues. PremiTest was able to detect all pen G residues in incurred samples containing pen G at concentrations greater than or equal to 12 ppb .

Table 9: Scanner results to confirm the Premi®Test's detection sensitivity as 10 ppb for penicillin G (1)

Zero control			
Sample	Interpreted Scanner Results	Sample	Interpreted Scanner Results
	NEG		
1	NEG	9	
2	NEG	10	
3	NEG	11	NEG
4	NEG	12	NEG
5	NEG	13	NEG
6	POS	14	NEG
7	NEG	15	NEG
8			NEG
Penicillin G 10 ppb			
Sample	Interpreted Scanner Results	Sample	Interpreted Scanner results
	POS		POS
	POS		POS
1	POS	16	POS
2	POS	17	POS
3	POS	18	POS
4	POS	19	POS
5	POS	20	POS
6	POS	21	POS
7	POS	22	POS
8	POS	23	POS
9	POS	24	POS
10	POS	25	POS
11	POS	26	POS
12	POS	27	POS
13	POS	28	POS
14		29	
15		30	

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Start	Complaint-No.	Lot	Complaint
26.03.2018	REM-R3900-115199	17H07/16	Liquid leaked
31.10.2018	REM-R3900-116382	18F01B16	Contents of the ampoules liquefied
02.01.2019	REM-R3900-116712	18F01/16	Deformations on the media, during pipetting
07.01.2019	REM-R3900-116732	18F01B16	Gel partially liquefied
25.01.2019	REM-R3925-116851	18F01/16	Tubes are insufficiently filled
10.04.2019	REM-R3925-117223	18L18/16	No color change
16.05.2019	REM-R3925-117374	18L10/16	Fill level
24.07.2019	REM-R3900-117745	18L18/16	2 Tubes are dried out
26.07.2019	REM-R3925-117755	19D15/16	Wrong delivery
16.08.2019	REM-R3925-117853	19A04/16	Delivery went to wrong customer
04.09.2019	REM-R3925-118114	19D15/16	Results for egg samples are not correct
15.10.2019	REM-R3900-118173	19B18/16	Results are not correct
18.10.2019	REM-R3925-118190	19D15/16	Ampoules underfilled
22.11.2019	REM-R3900-118372	19H12/16	Forwarder has not found address
20.07.2020	REM-R3900-119485	20E08/16	Customer complains about fill level
07.08.2020	REM-R3925-119550	20D15/16	underfilled ampoule
08.10.2020	REM-R3925-119842	20C10/16	displaced gel in ampoule
14.10.2020	REM-R3925-119863	20A20B16	Tube insufficiently filled
06.11.2020	REM-R3900-119954	20F11/16	Syringes for dosing not included in any kit
19/03/2021	REM-100277 - R3900	20I08/16	Delayed colour change
12/04/2021	REM-100367 - R3925	20I03/16	Gel partially liquefied
20/04/2021	REM-120059 - R3925	21A18/16	Component missing
28/05/2021	REM-120182 - R3900	21B25/16	no colour change for negative control after 3h
04/06/2021	REM-120198 - R3925	21B25/16	no colour change for negative control after 3h
30/06/2021	REM-120305 - R3925	21A18/16	Vial unfilled
27/07/2021	REM-120385 - R3925	21D14/16	discolored gel in a tube
25/08/2021	REM-120498 - R3925	21D13/16	Vial unfilled
16/12/2021	REM-120905 - R3900	21F16/16	no colour change for negative control after 3h
21/01/2022	REM-121002 - R3900	21F30/16	Component unfilled

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● Premi®Test detection limits in different animal food products

Substance	Chicken	Pork	Beef	Egg	Shrimp
β-lactams					
Amoxicillin	5	5	5	5	15
Ampicillin	5	5	5	5	
Penicillin-G	2,5	2,5	2,5	2,5	5
Cloxacillin		>100		100	
Oxacillin		100			
Dicloxacillin					
Cephalosporins					
Cefquinome	75	100	100		
Ceftiofur	100	200	100	400	
Macrolides					
Tylosin	50	25 - 50	50	50	
Erythromycin	100	100	100	50	100
Lincomycin	100	100	100		
Tilmicosin	50	50	50		
Spiramycin	1000	1000	1000		
Tetracyclines					
Chlortetracycline	100	100	100	600	1000
Oxytetracycline	100	100	100	400	100
Doxycycline	100	100	100	200	
Tetracycline		50		200	
Demeclocycline		50			
Sulphonamides					
Sulphamethazine	75	50-100	100	25	
Sulphadiazine	75	50-75	75	25	50
Sulphamethizole		50-100			
Sulphaguanidine	<200	150	<200		
Sulphadimethoxine		25 - 50	<100		50
Sulphapyridine	<50	50	<100		
Sulphamethoxypyridine	<100	25			
Sulphisoxazole	<100	25			
Sulphathiazole	<100	25			
Sulphachloropyridazine	<100	25			
Sulphamerazine	<100	25	<100		
Sulphanilamide	<100	150			
Sulphaquinoxaline	<100	50	<50		
Sulphamethiazole	<100		<50		
Sulfamethoxazole				25	

All detection limits are given in µg/kg = ppb

Substance	Chicken	Pork	Beef	Egg	Shrimp
Aminoglycosides					
Gentamicin	100	100	100	100	
Streptomycin	1500	1500	3000	1000	
Neomycine	300	300	300	300	200
Dihydrostreptomycin					
Quinolones					
Oxolinic acid					>10000
Enrofloxacin	>600	>600	>600		
Flumequine	>100	>100	>100		
Polypeptides					
Virginiamycin	500	500	500		
Bacitracin	500	500	500		
Zn-bacitracin	1250				
Colistin	>1000				
Ionophores					
Salinomycin	1000				
Monensin	1250				
Oligosaccharides					
Avilamycin	>5000				
Others					
Florphenicol	100	100	100		5000
Chloramphenicol	2500	2500	2500	2500	
Trimethoprim	50				
Narasin	1250				
Amprolium	>2000				
Phosphomycine	>1500				
Ronidazole					>5000
Furazolidone	>1500				