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- Dr. Sara Stead, Central Science Laboratory (CSL), Sand Hutton, York YO41 1LZ, England ;
- Dr. Steve Lehotay, USDA Agricultural Research Service (ARS), Eastern Regional Research Center, Microbial Biophysics and Residue Chemistry Research Unit, 600 East Mermaid Lane, Wyndmoor, PA 19038, USA.

Rapporteurs :

- Dr. Joe Boison, Canadian Food Inspection Agency, 116 Veterinary Road, S7N 2R3 Canada ;
- Dr. Michael Murphy, University of Minnesota.

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11. REFERENCES BIBLIOGRAPHIQUES

European Commission. 2002. Commission Decision (EC) N° 2002/657 of 12 August 2002 implementing Council Directive 96/23/EC concerning the performance of analytical methods and interpretation of results. Official Journal of European Communities. L221: 8-36.

Aboudi AR, Almashhadany DA, Hijazin M. 2022. Detection and impact of heat treatments on tetracycline residues in chicken meats sold in Amman city, Jordan. Food Research. 6:53-59.

Bahmani K, Shahbazi Y, Nikousefat Z. 2020. Monitoring and risk assessment of tetracycline residues in foods of animal origin. Food Science and Biotechnology. 29:441-448.

Berendsen BJA, Pikkemaat MG, Stolker LAM. 2011. Are antibiotic screening approaches sufficiently adequate? A proficiency test. Anal. Chim. Acta. 685:170-175.

Bor N, Seguino A, Sentamu DN, Chepyatich D, Akoko JM, Muinde P, Thomas LF. 2023. Prevalence of Antibiotic Residues in Pork in Kenya and the Potential of Using Gross Pathological Lesions as a Risk-Based Approach to Predict Residues in Meat. Antibiotics. 12:492.

Cantwell H, O'Keeffe M. 2006. Evaluation of the Premi Test and comparison with the One-Plate Test for the detection of antimicrobials in kidney. Food Additives and Contaminants. 23:120-125.

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Commission E. 2002. Commission Decision (EC) N° 2002/657 of 12 August 2002 implementing Council Directive 96/23/EC concerning the performance of analytical methods and interpretation of results. Off J Eur Comm. L221: 8-36.

CRL. 2010. Guideline for the validation of screening methods for residues of veterinary medicines (initial validation and transfer). Available from: https://food.ec.europa.eu/system/files/2016-10/cs_vet-med-residues_guideline_validation_screening_en.pdf. 1-18.

Dang PK, Degand G, Danyi S, Pierret G, Delahaut P, Ton VD, Maghuin-Rogister G, Scippo M-L. 2010. Validation of a two-plate microbiological method for screening antibiotic residues in shrimp tissue. Anal. Chim. Acta. 672:30-39.

Dang PK, Degand G, Douny C, Ton VD, Maghuin-Rogister G, Scippo ML. 2011. Optimisation of a new two-plate screening method for the detection of antibiotic residues in meat. Int. J. Food Sci. Technol. 46:2070-2076.

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Ekli R, Adzitey F, Huda N. 2019. Prevalence of resistant *Salmonella* spp. isolated from raw meat and liver of cattle in the Wa Municipality of Ghana. *IOP Conference Series: Earth and Environmental Science*. 287:012006.

Ezenduka EV, Ike OS, Anaelom NJ. 2014. Rapid detection of antimicrobial residues in poultry: A consequence of non-prudent use of antimicrobials. *Health*. 6:149-152.

Ezenduka EV, Okorie-Kanu OJ, Nwanta JA. 2019. Comparative analysis of two microbiological tests in the detection of oxytetracycline residue in chicken using ELISA as gold standard. *Journal of Immunoassay and Immunochemistry*. 40:617-629.

Fischer E, Huiñapi Z, Zegarra Paredes LF, Olortegui MP, Yori PP, Schiaffino F, Kosek MN. 2026. Prevalence of Antibiotic Residues in Commercial Poultry Feed and Meat from Industrial and Backyard Poultry in Iquitos, Peru. *IJID One Health*. 100112.

Gaudin V, Gauguain-Juhel M, Moretain JP, Sanders P. 2008. AFNOR validation of Premi®Test, a microbiological-based screening tube-test for the detection of antimicrobial residues in animal muscle tissue. *Food Addit. & Contam.: Part A*. 25:1451 - 1464.

Gaudin V, Hedou C, Rault A, Sanders P, Verdon E. 2009. Comparative study of three screening tests, two microbiological tube tests, and a multi-sulphonamide ELISA kit for the detection of antimicrobial and sulphonamide residues in eggs. *Food Addit. & Contam.: Part A*. 26:427-440.

Gondová Z, Kožárová I, Poláková Z, Mad'arová M. 2014. Comparison of Four Microbiological Inhibition Tests for the Screening of Antimicrobial Residues in the Tissues of Food-Producing Animals. *Italian Journal of Animal Science*. 13:3521.

Hoff R, Ribarcki F, Zancanaro I, Castellano L, Spier C, Barreto F, Fonseca SH. 2012. Bioactivity-based screening methods for antibiotics residues: A comparative study of commercial and inhouse developed kits. *Food Additives and Contaminants - Part A Chemistry, Analysis, Control, Exposure and Risk Assessment*. 29:577-586.

Kaba S, Koko KB, Zannou-Tchoko VJ, Konan BA, Datté JY. 2026. VETERINARY USE OF ANTIBIOTICS AND CONTAMINATION OF PORK CONSUMED IN CÔTE D'IVOIRE BY OXYTETRACYCLINE, TYLOSIN AND AMOXICILLIN RESIDUES. *Bulgarian Journal of Veterinary Medicine*. 29:180-188.

Kilinc B, Cakli S. 2008. Screening for antibiotic residues in the trout by the Four Plate test, Premi test and ELISA test. *Eur. Food Res. Technol*. 226:795-799.

Kilinc B, Meyer C, Hilge V. 2007. Evaluation of the EEC four-plate test and Premi test for screening antibiotic residues in trout (*Salmo trutta*). *Int. J. Food Sci. Technol*. 42:625-628.

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Kim DP, Degand G, Douny C, Pierret G, Delahaut P, Ton VD, Granier B, Scippo M-L. 2013. Preliminary Evaluation of Antimicrobial Residue Levels in Marketed Pork and Chicken Meat in the Red River Delta Region of Vietnam. Food and Public Health. 3 267-276.

Kirrella GAK, Deeb AMM, Abdallah RMI. 2017. Safety of frozen liver for human consumption. Journal of Food and Drug Analysis. 25:520-524.

Kozarova I, Janosova J, Mate D, Tkacikova S. 2009. Evaluation of three different microbial inhibition tests for the detection of sulphamethazine residues in the edible tissues of rabbit. Food Additives and Contaminants. Vol. 26,:978-987.

Kožárová I, Juščáková D, Šimková J, Milkovičová M, Kožár M. 2020. Effective screening of antibiotic and coccidiostat residues in food of animal origin by reliable broad-spectrum residue screening tests. Italian Journal of Animal Science. 19:487-501.

Kožárová I, Marcinčák S, Reitznerová A, Bartkovský M, Mačanga J, Marcinčáková D, Klempová T, Čertík M, *Verification for the presence of inhibitory substances in poultry meat after the consumption of the feed mixture supplemented with fermented feed.* 2016.

Magalhães CG, De Paiva CR, Botelho BG, De Oliveira AMG, De Souza LF, Nonaka CV, Santos KV, Farias LM, Carvalho MAR. 2012. In-house validation of Premi®Test, a microbiological screening test with solvent extraction, for the detection of antimicrobial residues in poultry muscles. Food Addit. & Contam.: Part A. 29:535-540.

Marcincak S, Hussein K, Mate D, Kozarova I, Sokol J, Zdolec N. 2006. Premi Test - a screening test for detecting sulphadimidine residues in eggs of laying hens. Medycyna Weterynaryjna. 62:159-161.

Mitiku B, Adugna M. 2018. Antimicrobial residue occurrence and its public health risk of beef meat in Debre Tabor and Bahir Dar, Northwest Ethiopia. Veterinary World. 11:902-908.

Mohamed MA, Ali A, Shair MA, Abatcha M, Hersi M, Ahmed A, Osman A. 2022. Detection of antibiotic residues in camel and goat meat from markets in Mogadishu city, Somalia. PAMJ - One Health. 9:

Mouliom Mouiche MM, Okah-Nnane NH, Moffo F, Djibo I, Mapiefou NP, Mpouam SE, Mfopit YM, Mingoas J-PK, Tebug SF, Ndukum JA. 2024. Antibiotic Residues in Foods of Animal Origin in Cameroon: Prevalence, Consumers' Risk Perceptions, and attitudes. J. Food Prot. 87:100237.

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Mustapha A, Musa S, Suleiman M, Abdullahi M, Isa T, Mahmoud F, Michael M, Mohammed J, Muhammad F, Balumi A, Ahmed U, Jacob L, Bukar H, Audu F, Musa M. 2023. Antibiotic Residues in fresh Chicken, Fish and Beef Meat samples in Maiduguri, Nigeria. Sahel Journal of Veterinary Sciences. 20:53-58.

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Nhung NT, Van NTB, Cuong NV, Duong TTQ, Nhat TT, Hang TTT, Nhi NTH, Kiet BT, Hien VB, Ngoc PT, Campbell J, Thwaites G, Carrique-Mas J. 2018. Antimicrobial residues and resistance against critically important antimicrobials in non-typhoidal *Salmonella* from meat sold at wet markets and supermarkets in Vietnam. *Int. J. Food Microbiol.* 266:301-309.

Njoga EO, Onunkwo JI, Okoli CE, Ugwuoke WI, Nwanta JA, Chah KF. 2018. Assessment of antimicrobial drug administration and antimicrobial residues in food animals in Enugu State, Nigeria. *Tropical Animal Health and Production.* 50:897-902.

Omeiza GK, Ajayi IE, Ode OJ, *Assessment of antimicrobial drug residues in beef in Abuja, the Federal Capital Territory, Nigeria.* 2012; Vol. 48, p 283-9.

Pikkemaat MG, Rapallini MLBA, Oostra-van Dijk S, Elferink JWA. 2009. Comparison of three microbial screening methods for antibiotics using routine monitoring samples. *Anal. Chim. Acta.* 637:298–304.

Pikkemaat MG, Rapallini MLBA, Zuidema T, Elferink JWA, Oostra-van Dijk S, Driessen-van Lankveld WDM. 2011. Screening methods for the detection of antibiotic residues in slaughter animals: comparison of the European Union Four-Plate Test, the Nouws Antibiotic Test and the PremiTest (applied to muscle and kidney). *Food Additives and Contaminants.* 28:26–34.

Rahimi Z, Shahbazi Y, Ahmadi F. 2018. Comparative Screening of Chloramphenicol Residue in Chicken Tissues Using Four Plate Test and Premi®Test Methods. *Pharm Sci.* 24:157-162.

Rahman MR, Lou Z, Antora R, Saqib M. 2015. Lateral Flow Test Strip Approaches For Rapid Detection Of Food Contaminants And Pathogens: Review. *American Research Thoughts.* 1:1222-1248.

Rakotoharinome M, Pognon D, Randriamparany T, Ming JC, Idoumbin J-P, Cardinale E, Porphyre V. 2014. Prevalence of antimicrobial residues in pork meat in Madagascar. *Tropical Animal Health and Production.* 46:49-55.

Sallam KI, Saad FSS, Abdelkhalek A. 2022. Health risk assessment of antimicrobial residues in sheep carcasses marketed in Kuwait. *Food Chem.* 383:132401.

Schneider M, Lehotay S. 2008. A comparison of the FAST, Premi® and KIS tests for screening antibiotic residues in beef kidney juice and serum. *Anal. Bioanal. Chem.* 390:1775-1779.

Schneider MJ, Mastovska K, Lehotay SJ, Lightfield AR, Kinsella B, Shultz CE. 2009. Comparison of screening methods for antibiotics in beef kidney juice and serum. *Anal. Chim. Acta.* 637:290–297.

Spišáková D, Kožárová I, Hriciková S, Marcinčák S. 2024. Comprehensive Screening of Salinomycin in Feed and Its Residues in Poultry Tissues Using Microbial Inhibition Tests Coupled to Enzyme-Linked Immunosorbent Assay (ELISA). *Foods.* 13:1661.

Staruch L, Mati M. 2012. Detection of inhibitory substances in poultry meat and poultry meat quality AND POULTRY MEAT QUALIT. *FOLIA VETERINARIA*, 56, Supplementum I: 51—53.

Tcheou PB, Bedekelabou AP, Adjei-Mensah B, Kpomasse CC, Talaki E, Salou M. 2025. Antibiotic use and residue detection in guinea fowl eggs in rural Togo: An assessment of practices and risks. *Veterinary and Animal Science*. 28:100448.

Vahedi Nouri N, Salehi A. 2024. Investigation of the antibiotic residues of broiler meat in northern Iran. *JSFA reports*. 4:33-38.

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Wu Q, Shabbir MAB, Peng D, Yuan Z, Wang Y. 2021. Microbiological inhibition-based method for screening and identifying of antibiotic residues in milk, chicken egg and honey. *Food Chem*. 363:130074.

Yang S-C, Yu M-C, Lee Y-H, Wang J-L, *Antibiotic Residues in Meat and Eggs in Taiwan: A Local Surveillance*. 2016; Vol. 12, p 1-6.

Ziani K, Pérez-López M, Mansouri A, Khaled MB, Rodriguez AS, Slimani M. 2018. Assessment of Oxytetracycline Residue in Cooked and Raw Meat of Chicken Broilers Before and After the End of Official Withdrawal Period. *Food Analytical Methods*. 11:2528-2537.

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No. Application: NF102

DECISION RELATED TO THE USAGE RIGHT

File followed by the Department: PGCP6
Certification Engineer: **Mélaine TERRO**
Phone: +33 (0)1 41 62 62 39
Email: melaine.terro@afnor.org

Operations supervisor : **Stella DUFLO**
Phone: 01 41 62 60 63
Email: stella.duflo@afnor.org

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La Plaine Saint-Denis, August 5th 2022

Dear Madam,

Following the positive agreement expressed on June 16th 2022, by the Technical Board of the NF VALIDATION mark (NF102), in its application to the food industry, I beg to inform you that the **NF VALIDATION certification has been renewed** for the following alternative method:

PREMI®TEST

Certificate reference No. RBP 31/02-04/11, with end of validity 30th-Aug.-2026

The alternative method has been validated for the detection of antibiotics residues in beef, pork and poultry meat (with the exception of ground meat), according to the validation protocol for detection and quantification of residues of veterinary drugs in food matrices - Rev. 1 (2017) of the NF VALIDATION mark, in its application to the food industry (NF102).

Yours Sincerely.

Mélaine TERRO
Certification Engineer

Annexe 2a. Notice d'utilisation 1011 (29/11/2010).

Premi®Test

RBP 31/02 - 04/11
ALTERNATIVE ANALYTICAL
METHODS FOR AGRIBUSINESS
Certified by AFNOR Certification
www.afnor-validation.org

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Premi®Test

Microbial Screening test for the presumptive detection of antibiotic residues in meat, fish, shrimps, eggs, kidneys, liver, urine, blood and feed. NFVALIDATION certification scope includes beef, pork and poultry meat (with the exception of ground meat).

RBP 31/02 - 04/11 ALTERNATIVE
ANALYTICAL METHODS FOR
AGRI-BUSINESS
Certified by AFNOR Certification
End of validity: 30th August 2018
<http://nf-validation.afnor.org>

2014/10/15 - insert 1011

Contents

25 ampoules with *Bacillus stearothermophilus* var. *calidolactis* in a solid agar medium, disposable pipette tips, perforated cover foil + syringe.

Product description

Premi®Test is a broad spectrum microbial screening test especially developed for the detection of antimicrobial substances, such as antibiotic and sulphonamide residues, in fresh meat, fish and eggs at or below most Maximum Residue Level (MRL).

Test principle

Premi®Test is based on the inhibition of the growth of *Bacillus stearothermophilus*, a micro-organism very sensitive to many antibiotics and sulphonamide residues.

A standardised number of spores are imbedded in an agar medium with selected nutrients. When meat fluid is added to the Premi®Test and heated at 64 °C, the spores will germinate. The germinated spores will multiply and form an acid when no inhibitory substances are present. This will be visible by a colour change from purple to yellow of the indicator in the ampoule. When antimicrobial residues are present in sufficient amount (above the detection level) no growth will occur and the colour will remain purple.

Special "sample-procedures" for fish, shrimps, eggs, kidneys, liver, urine, blood and feed are available at R-Biopharm AG.

Caution

This test is extremely sensitive to antibiotics and other inhibitory substances, any contamination with such materials should be avoided at all time. It is advised to wash hands thoroughly before starting the test procedure. Use tissue paper or a clean towel to dry hands.

Instruction for use

- Wash hands thoroughly before starting the test procedure.
- Cut out the required number of ampoules without damaging the aluminum foil of the adjacent ampoules.
- Remove foil carefully from the required amount of ampoules (do not open more ampoules as needed).
- Take approximately 2 cm³ of lean meat and use a meat press to extract approximately 250 µl of meat fluid. Alternatively, the meat juice can be obtained with the help of the Multipress or by freeze/thawing the meat. (A technical bulletin on the Multipress is available at www.r-biopharm.com).
- Use a new disposable tip on the syringe for each sample.
- Pipette 100 µl of the fluid onto the agar in the ampoule. Do not distort the agar.
- Allow to stand at room temperature for 20 minutes for a pre-incubation.
- Flush the meat juice away by gently filling and emptying the test ampoule twice with demineralized (demi) water. **Wash the ampoules with demineralized water only and do not use tap-water!**
- Remove the last water carefully from the test ampoule.
- Close the test ampoule with the foil supplied to avoid evaporation.
- Incubate the test ampoule in a Premi®Test Incubator or water bath (64 °C ± 0,5 °C).
- It is **strongly recommended** to use, in parallel, a negative control of the same matrix that is being tested. Read the test results after the negative control has changed color.
- All necessary equipment (meat press, scissors, incubator, and timer) are available in the Premi®Test Starter Kit.
- To facilitate the use of the Premi®Test there is a Powerpoint document available on request for all users.

Reading test results

- Read the color only from the bottom 2/3 part of the ampoule.
- A clear color change (from purple towards yellow) indicates the absence of antibiotics / sulphonamides above the detection limit.
- No clear color change indicates the presence of antibiotics and/or sulphonamides at/or above the limit of detection of the test.

Negative control

The use of a negative control is **strongly recommended (obligatory for NFVALIDATION)**. Start to check the color of the negative control after an incubation time of 2 h 40 min and then in intervals of 5 min until the color of the negative control has changed from purple to yellow. If there is no color change of the negative control after **3 h 10 min**, repeat the test. The negative control can be prepared by freezing and storing meat juice from a sample tested negative. **Never use water as a negative control!**

Positive control

It is recommended to test regularly a positive control (protocol is available from R-Biopharm AG) to verify the correct application of the test.

Storage

The ampoules should be stored cool (3 – 10 °C). Attention: **DO NOT FREEZE!**

Limited Liability

Premi®Test is a screening test and as such 100% accuracy of the test results cannot be guaranteed. Besides, the assessment of colour, in particular that of a yellow/purple result, may differ from person to person. In cases where severe consequences are involved for the user, test results should be confirmed by a validated comprehensive analytical method. The data corresponds to our present state of technology and provides information on our products and their uses. R-Biopharm makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. Defective products will be replaced. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. R-Biopharm shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

R-Biopharm AG
Postanschrift / Postal Address:
An der neuen Bergstraße 17
64297 Darmstadt, Germany
Sitz / Corporate Seat: Pfungstadt
Tel.: +49 (0) 61 51 - 81 02-0
Fax: +49 (0) 61 51 - 81 02-40
E-mail: Info@r-biopharm.de
www.r-biopharm.com

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Premi®Test

RBP 31/02-04/11
ALTERNATIVE ANALYTICAL
METHODS FOR AGRIBUSINESS
Certified by AFNOR Certification
End of validity: 30th August 2022
<http://nf-validation.afnor.org>

Microbial Screening test for the presumptive detection of antibiotic and sulfonamide residues in particularly fresh meat. NF VALIDATION certification scope includes beef, pork and poultry meat (with the exception of ground meat). The following selected antibiotics were investigated:
Penicillins: Penicillin G (CC β = 6 μ g/kg); Amoxicillin (CC β = 11 μ g/kg); Cloxacillin (CC β = 150 μ g/kg)
Tetracyclines: Chlortetracycline (CC β = 160 μ g/kg); Oxytetracycline (CC β = 160 μ g/kg); Doxycycline (CC β = 100 μ g/kg)
Sulfonamides: Sulfadimethoxine (CC β = 75 μ g/kg); Sulfadiazine (CC β = 90 μ g/kg)
Macrolides: Erythromycin A (CC β = 200 μ g/kg); Tylosine A (CC β = 90 μ g/kg)
The CC β values were determined by spiking in meat juice; the detection capabilities in other samples may differ.
Cephalosporins, Aminoglycosides, Amphenicols, Quinolones, Tiamulin, and Lincomycin are detected above the Maximum Residue Level (MRL).

2019/03/21 - Insert 1011

Contents

25 ampoules with *Bacillus stearothermophilus* var. *calidolactis* in a solid agar medium, disposable pipette tips, perforated cover foil + syringe.

Product description

Premi®Test is a broad spectrum microbial screening test especially developed for the detection of antimicrobial substances, such as antibiotic and sulfonamide residues particularly in fresh meat, in the range of MRL.

Test principle

Premi®Test is based on the inhibition of the growth of *Bacillus stearothermophilus*, a micro-organism very sensitive to many antibiotic and sulfonamide residues.

A standardised number of spores are imbedded in an agar medium with selected nutrients. When meat fluid is added to the Premi®Test and heated at 64 °C, the spores will germinate. The germinated spores will multiply and form an acid when no inhibitory substances are present. This will be visible by a colour change from purple to yellow of the indicator in the ampoule. When antimicrobial residues are present in sufficient amount (above the detection level) no growth will occur and the colour will remain purple.

Caution

This test is sensitive to certain antibiotics and other inhibitory substances, any contamination with such materials should be avoided at all time. It is advised to wash hands thoroughly before starting the test procedure. Use tissue paper or a clean towel to dry hands.

Instruction for use

- Wash hands thoroughly before starting the test procedure.
- Cut out the required number of ampoules without damaging the aluminum foil of the adjacent ampoules.
- Remove foil carefully from the required amount of ampoules (do not open more ampoules as needed).
- Take approximately 2 cm³ of lean meat and use a meat press to extract approximately 250 μ l of meat fluid. Alternatively, the meat juice can be obtained with the help of the Multipress or by freeze/thawing the meat. (A technical bulletin on the Multipress is available at www.r-biopharm.com).
- Use a new disposable tip on the syringe for each sample.
- Pipette 100 μ l of the fluid onto the agar in the ampoule. Do not distort the agar.
- Allow to stand at room temperature for 20 minutes for a pre-incubation.
- Flush the meat juice away by gently filling and emptying the test ampoule twice with demineralized (demi) water. Wash the ampoules with demineralized water only and do not use tap-water!
- Remove the last water carefully from the test ampoule.
- Close the test ampoule with the foil supplied to avoid evaporation.
- Incubate the test ampoule in a Premi®Test Incubator or water bath (64 °C $\pm$ 1 °C).
- It is required to use, in parallel, a negative control of the same matrix that is being tested. Read the test results after the negative control has changed color.
- All necessary equipment (meat press, scissors, incubator, and timer) are available in the Premi®Test Starter Kit.
- To facilitate the use of the Premi®Test there is a Powerpoint document available on request for all users.

Reading test results

- Read the color only from the bottom 2/3 part of the ampoule.
- A clear color change (from purple towards yellow) indicates the absence of antibiotics / sulfonamides above the detection limit.
- No clear color change indicates the presence of antibiotics and/or sulfonamides at/or above the limit of detection of the test.

Negative control

The use of a negative control is required (obligatory for NF VALIDATION). Start to check the color of the negative control after an incubation time of 2 h 40 min and then in intervals of 5 min until the color of the negative control has changed from purple to yellow. At that time, the samples should be compared to the negative control, because a delayed measuring (> 5 min.) can lead to increased values of the detection capability (CC β). If there is no color change of the negative control after 4 h, repeat the test. The negative control can be prepared by freezing and storing meat juice from a sample tested negative. Never use water as a negative control!

Positive control

It is strongly recommended to test a positive control (protocol is available from R-Biopharm AG) to verify the correct application of the test.

Storage

The ampoules should be stored cool (3 – 10 °C). Attention: DO NOT FREEZE!

Limited Liability

Premi®Test is a screening test and as such 100% accuracy of the test results cannot be guaranteed. Besides, the assessment of colour, in particular that of a yellow/purple result, may differ from person to person. In cases where severe consequences are involved for the user, test results should be confirmed by a validated comprehensive analytical method. R-Biopharm makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality.

Recommendations for sample-preparations of fish, shrimps, eggs, kidney, liver, urine, blood and feed are available at R-Biopharm AG. These matrices require individual validation by the customer. Supportive information are available at R-Biopharm upon request.

R-Biopharm AG
Postanschrift / Postal Address:
An der neuen Bergstraße 17
64297 Darmstadt, Germany
Sitz / Corporate Seat: Pfungstadt
Tel.: +49 (0) 61 51 - 81 02-0
Fax: +49 (0) 61 51 - 81 02-40
E-mail: info@r-biopharm.de
www.r-biopharm.com

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Premi®Test

RBP 31/02-04/11
ALTERNATIVE ANALYTICAL
METHODS FOR AGRIBUSINESS
Certified by AFNOR Certification
End of validity: see AFNOR website
<http://nf-validation.afnor.org/en>,
search PremiTest.

Microbial Screening test for the presumptive detection of antibiotic and sulfonamide residues in particularly fresh meat. NF VALIDATION certification scope includes beef, pork and poultry meat (with the exception of ground meat). The following selected antibiotics were investigated:
Penicillins: Penicillin G (CCB = 6 µg/kg); Amoxicillin (CCB = 11 µg/kg); Cloxacillin (CCB = 150 µg/kg)
Tetracyclines: Chlortetracycline (CCB = 160 µg/kg); Epichlortetracycline (CCB = 500 µg/kg); Oxytetracycline (CCB = 160 µg/kg); Epioxytetracycline (CCB = 500 µg/kg); Doxycycline (CCB = 100 µg/kg)
Sulfonamides: Sulfadimethoxine (CCB = 75 µg/kg); Sulfadiazine (CCB = 90 µg/kg)
Macrolides: Erythromycine A (CCB = 200 µg/kg); Tylosine A (CCB = 90 µg/kg)
The CCB values were determined by spiking in meat juice; the detection capabilities in other samples may differ.
Tetracyclines, Cefalosporins, Aminoglycosides, Amphenicols, Quinolones, Tiamulin, and Lincomycin are detected above the Maximum Residue Level (MRL).

20220516 - Insert 1012

Contents

25 ampoules with *Bacillus stearothermophilus* var. *calidolactis* in a solid agar medium, disposable pipette, perforated cover foil.

Product description

Premi®Test is a broad-spectrum microbial screening test especially developed for the detection of antimicrobial substances, such as antibiotic and sulfonamide residues particularly in fresh meat, in the range of MRL.

Test principle

Premi®Test is based on the inhibition of the growth of *Bacillus stearothermophilus*, a micro-organism very sensitive to many antibiotic and sulfonamide residues.

A standardized number of spores are imbedded in an agar medium with selected nutrients. When meat fluid is added to the Premi®Test and heated at 64 °C, the spores will germinate. The germinated spores will multiply and form an acid when no inhibitory substances are present. This will be visible by a color change from purple to yellow of the indicator in the ampoule. When antimicrobial residues are present in sufficient amount (above the detection level) no growth will occur and the color will remain purple.

Caution

This test is sensitive to certain antibiotics and other inhibitory substances, any contamination with such materials should be avoided at all times. It is advised to wash hands thoroughly before starting the test procedure. Use tissue paper or a clean towel to dry hands.

Instruction for use

- Wash hands thoroughly before starting the test procedure.
- Remove the required number of ampoules without damaging the aluminum foil of the adjacent ampoules.
- Remove foil carefully from the required number of ampoules (do not open more ampoules than needed).
- Take approximately 2 cm³ of lean meat and use a meat press to extract approximately 250 µL of meat fluid. Alternatively, the meat juice can be obtained by freeze/thawing the meat.
- For each sample use a new and clean pipette.
- Pipette 100 µL of the fluid onto the agar in the ampoule. Do not distort the agar.
- Allow to stand at room temperature for 20 minutes for a pre-incubation.
- Flush the meat juice away by gently filling and emptying the test ampoule twice with demineralized (demi) water. Wash the ampoules with demineralized water only and do not use tap-water!
- Remove the last water carefully from the test ampoule.
- Close the test ampoule with the foil supplied to avoid evaporation.
- Incubate the test ampoule in a Premi®Test Incubator or water bath (64 °C ± 1 °C).
- It is required to use, in parallel, a negative control of the same matrix that is being tested. Read the test results after the negative control has changed color.
- All necessary equipment (meat press, scissors, incubator, and timer) are available in the Premi®Test Starter Kit.
- To facilitate the use of the Premi®Test there is a PowerPoint document available on request for all users.

Reading test results

- Read the color only from the bottom 2/3 part of the ampoule.
- A clear color change (from purple towards yellow) indicates the absence of antibiotics / sulfonamides above the detection limit.
- No clear color change indicates the presence of antibiotics and/or sulfonamides at/or above the limit of detection of the test.

Negative control

The use of a negative control is required (obligatory for NF VALIDATION). Start to check the color of the negative control after an incubation time of 2 h 40 min and then in intervals of 5 min until the color of the negative control has changed from purple to yellow. At that time, the samples should be compared to the negative control, because a delayed measuring (5 min.) can lead to increased values of the detection capability (CCB). If there is no color change of the negative control after 4 h, repeat the test. The negative control can be prepared by freezing and storing meat juice from a sample tested negative. Never use water as a negative control!

Positive control

It is strongly recommended to test a positive control (protocol is available from R-Biopharm AG) to verify the correct application of the test.

Storage

The ampoules should be stored cool (3 – 10 °C). Attention: DO NOT FREEZE!

Limited Liability

Premi®Test is a screening test and as such 100% accuracy of the test results cannot be guaranteed. Besides, the assessment of color, in particular that of a yellow/purple result, may differ from person to person. In cases where severe consequences are involved for the user, test results should be confirmed by a validated comprehensive analytical method. R-Biopharm makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality.

Recommendations for sample-preparations of fish, shrimps, eggs, kidney, liver, urine, blood and feed are available at R-Biopharm AG. These matrices require individual validation by the customer. Supportive information are available at R-Biopharm upon request.

R-Biopharm AG
Postanschrift / Postal Address:
An der neuen Bergstraße 17
64297 Darmstadt, Germany
Sitz / Corporate Seat: Pfungstadt
Tel.: +49 (0) 61 51 - 81 02-0
Fax: +49 (0) 61 51 - 81 02-40
E-mail: info@r-biopharm.de
www.r-biopharm.com

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Premi®Test

RBP 31/02-04/11
ALTERNATIVE ANALYTICAL
METHODS FOR AGRIBUSINESS
Certified by AFNOR Certification
Fin de validité: see AFNOR website
<http://nf-validation.afnor.org/en>,
search PremiTest.

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Microbial Screening test for the presumptive detection of antibiotic and sulfonamide residues in particularly fresh meat. NFVALIDATION certification scope includes beef, pork and poultry meat (with the exception of ground meat). The following selected antibiotics were investigated:
Penicillins: Penicillin G (CCB = 6 µg/kg); Amoxicillin (CCB = 11 µg/kg); Cloxacillin (CCB = 150 µg/kg)
Tetracyclines: Chlortetracycline (CCB = 160 µg/kg); Epichlortetracycline (CCB = 500 µg/kg); Oxytetracycline (CCB = 160 µg/kg); Epioxytetracycline (CCB = 500 µg/kg); Doxycycline (CCB = 100 µg/kg)
Sulfonamides: Sulfadimethoxine (CCB = 75 µg/kg); Sulfadiazine (CCB = 90 µg/kg)
Macrolides: Erythromycin A (CCB = 200 µg/kg); Tylosine A (CCB = 90 µg/kg)
The CCB values were determined by spiking in meat juice; the detection capabilities in other samples may differ.
Tetracyclines, Cephalosporins, Aminoglycosides, Amphenicols, Quinolones, Tiamulin, and Lincomycin are detected above the Maximum Residue Level (MRL).

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Reading test results

- Read the color only from the bottom 2/3 part of the ampoule.
- A clear color change (from purple towards yellow) indicates the absence of antibiotics / sulfonamides above the detection limit.
- No clear color change indicates the presence of antibiotics and/or sulfonamides at/ or above the limit of detection of the test.

Negative control

The use of a negative control is required (obligatory for NFVALIDATION). Start to check the color of the negative control after an incubation time of 2 h 40 min and then in intervals of 5 min until the color of the negative control has changed from purple to yellow. At that time, the samples should be compared to the negative control, because a delayed measuring (> 5 min.) can lead to increased values of the detection capability (CCB). If there is no color change of the negative control after 4 h, repeat the test. The negative control can be prepared by freezing and storing meat juice from a sample tested negative. Never use water as a negative control!

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Annexe 2i. Notice d'utilisation 1014 (21/05/2026).

Premi®Test

ISO 16343-2:2015
ALTERNATIVE ANALYTICAL
METHODS FOR AGRI-FOODS
Certified by IFRC Certification
Not of validity: see IFRC's website
<https://ifrc-certification.org/en/our-services>
search Premi®Test

Microbial Screening test for the presumptive detection of antimicrobial residues in particularly fresh meat. NF VALIDATION certification scope includes beef, pork and poultry meat (after the evisceration of ground meat). The following selected antibiotics were investigated:
Penicillins: Penicillin G (CCB = 6 µg/kg), Amoxicillin (CCB = 11 µg/kg), Cloxacillin (CCB = 183 µg/kg)
Tetracyclines: Chlorotetracycline (CCB = 100 µg/kg), Tetracycline (CCB = 100 µg/kg), Doxycycline (CCB = 183 µg/kg), Erythromycin (CCB = 183 µg/kg)
Sulfonamides: Sulfadiazine (CCB = 75 µg/kg), Sulfamethoxazole (CCB = 90 µg/kg)
Macrolides: Erythromycin A (CCB = 100 µg/kg), Tylosine A (CCB = 90 µg/kg)
The CCB values were determined by spiking in meat (value, the detection capabilities in other samples may differ).
Tetracyclines, Cephalosporins, Aminoglycosides, Amphenicols, Quinolones, Fluoroquinolones, and Lincosamides are detected above the Maximum Residue Limit (MRL).
This validation study has been conducted according to the protocol NF ISO 16343-2:2015 for validation protocol for methods of detecting and quantifying traces of animal residues in food products (Ann. 12).

20200221 / Issue 03/IV

Conditions

25 ampoules with 100 µl of water (Premi®Test 100 µl) in a solid agar medium, **agarose plate**, performed over 6h.

Product description

Premi®Test is a broad spectrum microbial screening test especially developed for the detection of antimicrobial substances, such as antibiotics and sulfonamide residues particularly in fresh meat, in the range of MRL.

Test principle

Premi®Test is based on the inhibition of the growth of *Escherichia coli* (ATCC 8739), a micro-organism very sensitive to many antibiotics and sulfonamide residues.
A standardized number of spores are inoculated in an agar medium with selected nutrients. When meat fluid is added to the Premi®Test and heated at 65 °C, the spores will germinate. The germinated spores will multiply and form an acid when no inhibitory substances are present. This will be visible by a color change from purple to yellow of the indicator in the ampoule. When antimicrobial residues are present in sufficient amount (above the detection limit) no growth will occur and the color will remain purple.

Caution

This test is sensitive to certain antibiotics and other inhibitors substances, any contamination with such materials should be avoided at all times. It is advised to wash hands thoroughly before starting the test procedure. Use tissue paper or a clean towel to dry hands.

Instructions for use

- Wash hands thoroughly before starting the test procedure.
- Remove the required number of ampoules without damaging the aluminum foil of the adjacent ampoules.
- Remove foil carefully from the required number of ampoules (do not open more ampoules than needed).
- Take approximately 2 cm of meat (not a meat piece) and cut it into small pieces (approx. 2x2 cm) (alternatively, the meat juice can be obtained by pressing/tearing the meat).
- For each sample use a new and clean pipette.
- Pipette 100 µl of the fluid onto the agar in the ampoule. Do not disturb the agar.
- Allow to stand at room temperature for 20 minutes for a pre incubation.
- Push the meat juice evenly by gently filling and swirling the test ampoule twice with deiontized (dist.) water. Wash the ampoules with deiontized water only and do not use tap water!
- Remove the test water carefully from the test ampoule.
- Close the test ampoule with the foil applied to avoid evaporation.
- Incubate the test ampoule in a Premi®Test incubator or water bath (65 °C ± 1 °C).
- It is essential to use, in parallel, a negative control of the same matrix that is being tested. Read the test results after the negative control has changed color.
- All necessary equipment (meat press, scissors, incubator, and timer) are available in the Premi®Test Starter Kit.
- To facilitate the use of the Premi®Test there is a PaperPoint document available on request for all users.

Reading test results

- Read the color only from the bottom 1/3 part of the ampoule.
- After color change from purple towards yellow indicates the absence of antibiotic / sulfonamides above the detection limit.
- No color change indicates the presence of antibiotic and/or sulfonamides and/or above the limit of detection of the test.

Negative control

The use of a negative control is required (obligatory for NF VALIDATION). Start to check the color of the negative control after an incubation time of 2 h 40 min and then in intervals of 1 min until the color of the negative control has changed from purple to yellow. At that time, the samples should be compared to the negative control, because a delayed measuring (≥ 6 min) can lead to increased values of the detection capability (CCB). If there is no color change of the negative control after 4h, repeat the test. The negative control can be prepared by freezing and thawing meat juice from a sample tested negative. Never use water as a negative control.

Positive control

It is strongly recommended to test a positive control (provided) is available from R-Biopharm AG) to verify the correct application of the test.

Storage

The ampoules should be stored cool (5 - 10 °C). Attention: DO NOT FREEZE!


Limited Liability

Premi®Test is a screening test and as such 100% accuracy of the test results cannot be guaranteed. Besides, the assessment of color, in particular that of a yellow/purple result, may differ from person to person. In cases where severe consequences are involved for the user, test results should be confirmed by a validated comprehensive analytical method. R-Biopharm makes no warranty of any kind, whether expressed or implied, except that the materials from which its products are made are of standard quality.

Recommendations for sample preparations of fish, shrimp, eggs, honey, beer, urine, blood and feed are available at R-Biopharm AG. These materials require individual validation by the customer. Supportive information are available at R-Biopharm upon request.

R-Biopharm AG
Postfach 101 / Postal Address: An
der Leber 101
63071 Darmstadt, Germany
800 / Corporate line: Pfungstadt
Tel: +49 (0) 61 81 - 81 22 13
Fax: +49 (0) 61 81 - 81 22 40
E-mail: info@rbiopharm.com
www.rbiopharm.com

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Premi®Test

Microbial inhibition test validated for the detection of residues of PENICILLIN in fresh bovine meat

2022/01/7 issue/0007

1. The Intended Use

Premi®Test is a Microbial inhibition test validated for the detection of residues of PENICILLIN in fresh bovine meat.

2. The Intended User

The test kit is designed for use by quality control personnel and others familiar with the antibiotic screening testing. All users of the Premi®Test should be trained by a Biopharm AG representative or by somebody who already successfully completed the training. To facilitate the use of the Premi®Test there is a French/English document available for all users.

3. Contents

- 20 ampoules containing an agar medium, pH indicator and spores of *Bacillus Stenothermophilus*
- 20 single use pipettes
- Premi®Test covering kit

4. Materials required but not included

- Premi®Test Starter Kit, containing the required equipment to do Premi®Test
- Biopharm AG heating block (Premi®Test incubator), meat press, pair of scissors, thermometer and timer.
- Instead of the Premi®Test incubator a water bath at 65 °C ± 1.0 °C can be used.
- Instead a water bath is used a rack or flask is required to accommodate the ampoules during incubation.
- Calibrated thermometer appropriate for designated use.
- Refrigerator at 0 - 10 °C (32 - 50 °F) for storage of the Premi®Test ampoules.
- Negative control muscle tissue (juice from a known antibiotic free sample, previously tested as negative).
- Premi®Scan, software to interpret the results of the Premi®Test, an HP (Gigaset 7620 or 8230 scanner), template, black cloth and Excel 2002 Color chart.

5. Test Performance

5.1 Detection sensitivity (false response) (Premi®Test)

32 ppb for Penicillin G

The Premi®Test has been validated to detect penicillin G residues in fresh bovine muscle tissue consistently with 99% confidence.

5.2 Selectivity

32 negative control meat samples were evaluated in an independent laboratory and 3 of these negative samples tested positive with the Premi®Test.

5.2.1 Cross reactivity

Premi®Test does cross react with the following drugs spiked at a concentration of 100 ppb

β-lactams: Amoxicillin, Ampicillin, Cephalixin, Cloxacillin

Tetracyclines: Chloramphenicol, Tetracycline

Macrolides: Erythromycin, Tylosin, Thiomycin

Glycopeptides: Vancomycin

Sulfonamides: Sulfadiazine, Sulfadimethoxine, Sulfathiazole

Misfoldamycin: Capreomycin

Premi®Test does cross react with the following drugs spiked at a concentration of 50 ppb

β-lactams: Amoxicillin, Ampicillin, Cephalixin

Macrolides: Erythromycin, Tylosin, Thiomycin

Misfoldamycin: Capreomycin

Premi®Test does cross react with the following drugs spiked at a concentration of 10 ppb

β-lactams: Amoxicillin, Ampicillin

Misfoldamycin: Capreomycin

6. Test principle

- Premi®Test is based on the inhibition of the growth of *Bacillus Stenothermophilus* var. *colliditellus*, a thermophilic microorganism which is very sensitive to many antibiotics and sulfonamide residues. A standardized number of spores is inoculated in an agar medium with selected nutrients and a pH indicator (from neutral purple).
- When muscle tissue fluid containing no antimicrobials or non-inhibitory substances is added to the Premi®Test and heated at 65 °C, the spores will grow, multiply and produce acid. This process is characterized by a visible change in the color of the pH indicator from purple to yellow in the ampoule.
- When muscle tissue fluid containing antimicrobial residues or inhibitory substances are present at or above the detection sensitivity, no growth will occur and there will be no acid production. As a result, there will be no change in the color of the pH indicator and the color of the ampoule will remain purple.

7. Cautions

This test is extremely sensitive to antibiotics, disinfectants and other inhibitory substances. Any contamination with such materials should be avoided at all times. It is advised to wash hands thoroughly before starting the test procedure, and to dry the hands with tissue paper or a clean towel. Disposable hygienic hand gloves can also be used.

8. Operating instructions Premi®Test

8.1 Before starting the test procedure

1. Premi®Test: take the number of ampoules needed from storage and let them warm up to room temperature (20 - 25 °C / 64 - 77 °F)

2. Premi®Test incubator: Allow the Premi®Test incubator or water bath to reach the required temperature of 65 °C ± 1.0 °C (147 °F).

3. Samples: Take approximately 4 grams of fresh, lean bovine muscle tissue and use a meat press to extract approximately 250 % of meat fluid. The juice can be obtained somewhat easier, when heating the meat press first (e.g. by placing it in an incubator of 55 °C (132 °F)).

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4. Number of samples: The number of meat samples that can be tested in one series is 5. The Premi@Test incubator contains 10 holes. One hole should be reserved for the negative control.
5. Wash hands thoroughly before starting the test procedure.
- 5.2 Testing procedure**
 1. Remove the required number of ampoules. Do not tear off ampoules, as this may damage the foil of adjacent ampoules. Any test with torn foil must be used immediately or discarded.
 2. Mark the ampoules with a meat sample identification code.
 3. Remove foil carefully from the required amount of ampoules (do not open more ampoules than needed after opening, the test can not be stored any longer).
 4. Pick up a disposable pipette tip without touching the ends and place it on the syringe.
 5. For each sample use a new and clean pipette.
 6. Pipette 100 µl of 1% of each muscle tissue juice sample to be tested into the labeled ampoules. Depress the plunger of the pipette, dip the pipette tip into muscle tissue juice and let the plunger return slowly to its original position. Transfer the meat juice onto the agar of the labeled ampoules by depressing the plunger slowly and completely. Do not disturb the agar.
 7. Pre-incubate the test 20 minutes at room temperature (18 - 25 °C / 64 - 77 °F).
 8. Remove the muscle tissue juice by gently emptying the test, by putting the ampoule(s) upside down.
 9. Wash the test base lines, by gently filling and emptying the test ampoules with deionized (distilled) water.
 10. Remove the test bases of water as fully from the test ampoules, by putting the ampoules upside down on a tissue.
 11. Close the test ampoules with the foil supplied with the kit to avoid evaporation during incubation.
 12. Incubate the test ampoules in the Premi@Test incubator at water bath (55 °C / 131 °F), until the negative control changes color from purple to yellow (See Section 5.3).
 13. When the negative control has changed color (after approximately 5 hours), the ampoules can be read with the Premi@Scan.
- 5.3 Guidelines to Determine the Readout of Samples from the Incubator**
 1. Interpreting results for screening to meat color changed. **NOTE:** The VALIDATION WAS NOT INTENDED TO COVER THE USE OF THE PREM@TEST as a guideline for triggering the point at which the test samples should be removed from the incubator to be read by the Premi@Scan. **DISCLAIMER FOR TRANSMISSION OF THE POINT AT WHICH THE TEST SAMPLES SHOULD BE REMOVED FROM THE INCUBATOR TO BE READ BY THE PREMI@SCAN.**
 2. Read the color only from the bottom 20 mm of the ampoules.
 3. A partly or completely yellow color of the agar medium indicates the absence of antibiotic residues in the corresponding meat sample as it shows the weak indicated above. Such a result should be regarded as a negative result (important as not found).
 4. When the color changes to "doublet", incubate 5 minutes longer and read again. When the result is still not clear, test the sample again.
- 5. Operating Instructions Premi@Scan**
 1. Place the Premi@Scan in the CD drive; the program will start up automatically.
 2. Accept the license agreement by clicking on OK and the program will be ready to use.
 3. Place the Serial DSC Color chart, delivered with the Premi@Scan, on the scanner (50 Samples 7400 or 8200).
 4. Calibrate the Premi@Scan with the Serial DSC Color chart, by following the instructions of the software.
 5. Put the ampoules to be tested in the templates, delivered with the Premi@Scan.
 6. Place the black cloth over the templates.
 7. The Premi@Scan recognizes the position of the ampoules automatically and determines the result of the ampoules.
 8. Results can be printed.
- 5B Storage and Precautions**
 1. The Premi@Test ampoules should be stored in a dry and cool place (preferable below 5 °C / 41 °F). Attention: DO NOT FREEZE ampoules should not be used after the expiration date, indicated on the label.
 2. Anything that may come in contact with the meat (hands, pipette tips, etc.) must be cleaned thoroughly, rinsed with tap water and dried with a clean (paper) towel to avoid contamination with inhibiting substances like antibiotics, other drugs, cleaning agents or disinfectants.
 3. Store the meat below 5 °C (41 °F). If it is not possible to carry out the test immediately, it is not possible to test the meat within 48 hours, the sample should be stored at or below -18 °C.
- 16 Technical Assistance**

For assistance in use, interpretation, or troubleshooting contact the results of this test contact sales@r-biopharm.de. To facilitate the use of Premi@Test there is a Premi@Point document available for all users.
- 18 Limited Liability**


Premi@Test is a screening test and as such 100% accuracy of the test results cannot be guaranteed. Besides, the assessment of color, in particular that of a yellow/brown result, may differ from person to person. In cases where severe consequences are involved for the user, test results should be confirmed by a validated comprehensive analytical method (like HPLC or GC-MS/MS). The data corresponds to our present state of technology and provides information on our products and their uses. R-Biopharm makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. Defective products will be replaced. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. R-Biopharm shall not be liable for any damages, including special or consequential damages, or expenses arising directly or indirectly from the use of this product. The above text has received AGAC approval for this test kit use in North America.

R-Biopharm AG
Postfach 18 (Postal Address: An
der Leber 18) 717
65021 Darmstadt, Germany
RBL Corporate Inc. (Pfungstadt)
Tel.: +49 (0) 61 81 81 81 82 40
Fax: +49 (0) 61 81 81 82 40
E-mail: sales@r-biopharm.de
www.r-biopharm.de

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Annexe 3a. Certificat No. 060601 d'Octobre 2021.

CERTIFICATION
AOAC® Performance TestedSM
Certificate No. 060601
The AOAC Research Institute hereby certifies the performance of the method known as:
Premi®Test
manufactured by DSM PREMI® Test B. V. P. O. Box 6500 6401 JH HEERLEN The Netherlands
This method has been evaluated in the AOAC® Performance Tested MethodsSM Program and found to perform as stated by the manufacturer contingent to the comments contained in the manuscript. This certificate means that an AOAC® Certification Mark License Agreement has been executed which authorizes the manufacturer to display the AOAC Performance TestedSM certification mark along with the statement - "THIS METHOD'S PERFORMANCE WAS REVIEWED BY AOAC RESEARCH INSTITUTE AND WAS FOUND TO PERFORM TO THE MANUFACTURER'S SPECIFICATIONS" - on the above-mentioned method for a period of one calendar year from the date of this certificate (October 23, 2021 – December 31, 2022). Renewal may be granted at the end of one year under the rules stated in the licensing agreement.
 _____ Scott Coates, Senior Director Signature for AOAC Research Institute October 23, 2021 _____ Date
2275 Research Blvd., Ste. 300, Rockville, Maryland, USA Telephone: +1-301-924-7077 Fax: +1-301-924-7089 Internet e-mail: aoacri@aoac.org * World Wide Web Site: http://www.aoac.org

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METHOD AUTHORS Angelique de Rijk	MANUFACTURER DSM Nutritional Products DSM Premi® Test BV P. O. Box 6500 6401 JH Heerlen The Netherlands	CERTIFICATION MARK LICENSEE R-Biopharm AG An der Neuen Bergstraße 17, 64297 Darmstadt, Germany
METHOD NAME Premi® Test	CATALOG NUMBERS R3925, R3900	
INDEPENDENT LABORATORY Central Science Laboratory (CSL) York YO41 1LZ, England USDA Agricultural Research Service Eastern Regional Research Center Microbial Biophysics and Residue Chemistry Research Unit 600 East Mermaid Lane Wyndmoor, PA 19038, USA	AOAC EXPERTS AND PEER REVIEWERS Joe Boison ¹ , Michael Murphy ² ¹ Canadian Food Inspection Agency, Saskatoon, SK, Canada ² University of Minnesota, Minneapolis, MN, USA	
APPLICABILITY OF METHOD Target analyte – Penicillin G Matrix – Bovine muscle tissue Performance claims – See Table 1		
ORIGINAL CERTIFICATION DATE September 27, 2006	CERTIFICATION RENEWAL RECORD Renewed annually through December 2022.	
METHOD MODIFICATION RECORD 1. November 2019 Level 1	SUMMARY OF MODIFICATION 1. Editorial changes to inserts.	
Under this AOAC® Performance Tested SM License Number, 060601 this method is distributed by: NONE	Under this AOAC® Performance Tested SM License Number, 060601 this method is distributed as: NONE	

PRINCIPLE OF THE METHOD (1)

The Premi®test is based on the inhibition of the growth of *Bacillus Stearothermophilus*, a bacterium which is very sensitive to many antibiotics and sulfa compounds. A standardized number of spores are imbedded in an agar medium with selected nutrients. A small amount of meat juice is added to the test tube. After a pre-incubation period of 20 minutes at room temperature (18-22 °C i.e. 64.4-71.6 °F), the test tube is washed and incubated for approximately 3 hours at 64°C (147°F). At this temperature, the spores will start to grow (germinate) and multiply with the production of acid if no antimicrobial (non-inhibitory) substances are present in the meat juice. This will be visible by a color change in the test tube from purple to yellow. When antimicrobial (inhibitory) substances are present (at or above the detection sensitivity), no growth will occur and the color in the test tube will remain purple.

DISCUSSION OF THE VALIDATION STUDY (1)

The results show that PremiTest met all the requirements for this AOAC validation study for the detection of penicillin G residues in bovine muscle tissue. The results also show that there is cross reactivity with many other antimicrobial compounds. PremiTest is a broad-spectrum screening test which can detect many relevant antimicrobials in meat.

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
	
CERTIFICATION	
AOAC Research Institute <i>Performance Tested MethodsSM</i>	
Certificate No. 060601	
The AOAC Research Institute hereby certifies the method known as:	
Premi®Test manufactured by DSM PREMI® Test B. V. P. O. Box 6500 6401 JH HEERLEN The Netherlands	
This method has been evaluated in the AOAC Research Institute <i>Performance Tested MethodsSM</i> Program and found to perform as stated in the applicability of the method. This certificate indicates an AOAC Research Institute Certification Mark License Agreement has been executed which authorizes the manufacturer to display the AOAC Research Institute <i>Performance Tested MethodsSM</i> certification mark on the above-mentioned method for the period below. Renewal may be granted by the Expiration Date under the rules stated in the licensing agreement.	
 _____ Scott Coates, Senior Director Signature for AOAC Research Institute	Issue Date December 19, 2022 Expiration Date December 31, 2023
2275 Research Blvd., Ste. 300, Rockville, Maryland, USA Telephone: +1-301-924-7077 Fax: +1-301-924-7089 Internet e-mail: aoacri@aoac.org * World Wide Web Site: http://www.aoac.org	

CERTIFICATION


AOAC Research Institute
Performance Tested MethodsSM

Certificate No.
060601

The AOAC Research Institute hereby certifies the method known as:

Premi®Test
manufactured by
DSM PREMI® Test B. V.
P. O. Box 6500
6401 JH HEERLEN
The Netherlands

This method has been evaluated in the AOAC Research Institute *Performance Tested MethodsSM* Program and found to perform as stated in the applicability of the method. This certificate indicates an AOAC Research Institute Certification Mark License Agreement has been executed which authorizes the manufacturer to display the AOAC Research Institute *Performance Tested MethodsSM* certification mark on the above-mentioned method for the period below. Renewal may be granted by the Expiration Date under the rules stated in the licensing agreement.

Bradley A. Stawick, Senior Director
Signature for AOAC Research Institute

Issue Date	January 17, 2024
Expiration Date	December 31, 2024

2275 Research Blvd., Ste. 300, Rockville, Maryland, USA Telephone: +1-301-924-7077 Fax: +1-301-924-7089
Internet e-mail: aoacri@aoac.org * World Wide Web Site: <http://www.aoac.org>

		
CERTIFICATION		
AOAC Research Institute <i>Performance Tested MethodsSM</i>		
Certificate No. 060601		
The AOAC Research Institute hereby certifies the method known as:		
Premi®Test manufactured for R-Biopharm AG An der neuen Bergstr. 17 64297 Darmstadt, Germany		
This method has been evaluated and certified according to the policies and procedures of the AOAC <i>Performance Tested MethodsSM</i> Program. This certificate indicates an AOAC Research Institute Certification Mark License Agreement has been executed which authorizes the manufacturer to display the AOAC Research Institute <i>Performance Tested MethodsSM</i> certification mark on the above-mentioned method for the period below. Renewal may be granted by the Expiration Date under the rules stated in the licensing agreement.		
 Bradley A. Stawick, Senior Director Signature for AOAC Research Institute	Issue Date Expiration Date	December 13, 2024 December 31, 2025
2275 Research Blvd., Ste. 300, Rockville, Maryland, USA Telephone: +1-301-924-7077 www.aoac.org		

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	
CERTIFICATION	
AOAC Research Institute <i>Performance Tested MethodsSM</i>	
Certificate No. 060601	
The AOAC Research Institute hereby certifies the method known as PremiTest manufactured by R-Biopharm AG an der neuen Bergstraße 64297 Darmstadt Germany	
This method has been evaluated and certified according to the policies and procedures of the AOAC <i>Performance Tested MethodsSM</i> Program. This certificate indicates an AOAC Research Institute Certification Mark License Agreement has been executed which authorizes the manufacturer to display the AOAC Research Institute <i>Performance Tested MethodsSM</i> certification mark on the above-mentioned method for the period below. Renewal may be granted by the Expiration Date under the rules stated in the licensing agreement.	
 _____ Bradley A. Stawick, AOAC Research Institute Senior Director	Issue Date January 14, 2026 Expiration Date December 31, 2026
2275 Research Blvd., Ste. 300, Rockville, Maryland, USA Telephone: +1-301-924-7077 www.aoc.org	

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AUTHORS Angelique de Rijk	MANUFACTURED FOR R-Biopharm AG An der Neuen Bergstraße 17, 64297 Darmstadt, Germany
METHOD NAME Premi® Test	CATALOG NUMBERS R3925, R3900
INDEPENDENT LABORATORY Central Science Laboratory (CSL) York YO41 1LZ, England USDA Agricultural Research Service Eastern Regional Research Center Microbial Biophysics and Residue Chemistry Research Unit 600 East Mermaid Lane Wyndmoor, PA 19038, USA	APPLICABILITY OF METHOD Target analyte – Penicillin G Matrix – Bovine muscle tissue Performance claims – See Table 1
ORIGINAL CERTIFICATION DATE September 27, 2006	CERTIFICATION RENEWAL RECORD Renewed through December 2026.
METHOD MODIFICATION RECORD 1. November 2019 Level 1 2. December 2022 Level 1 3. January 2024 Level 1 4. December 2024 Level 1	SUMMARY OF MODIFICATION 1. Editorial changes. 2. Editorial changes. 3. Editorial changes. 4. Editorial changes.
Under this AOAC <i>Performance Tested Methods</i> SM License Number, 060601 this method is distributed by: NONE	Under this AOAC <i>Performance Tested Methods</i> SM License Number, 060601 this method is distributed as: NONE

PRINCIPLE OF THE METHOD (1)

The Premi®test is based on the inhibition of the growth of *Bacillus Stearothermophilus*, a bacterium which is very sensitive to many antibiotics and sulfa compounds. A standardized number of spores are imbedded in an agar medium with selected nutrients.

A small amount of meat juice is added to the test tube. After a pre-incubation period of 20 minutes at room temperature (18-22 °C i.e. 64.4-71.6 °F), the test tube is washed and incubated for approximately 3 hours at 64°C (147°F).

At this temperature, the spores will start to grow (germinate) and multiply with the production of acid if no antimicrobial (non-inhibitory) substances are present in the meat juice. This will be visible by a color change in the test tube from purple to yellow. When antimicrobial (inhibitory) substances are present (at or above the detection sensitivity), no growth will occur and the color in the test tube will remain purple.

DISCUSSION OF THE VALIDATION STUDY (1)

The results show that PremiTest met all the requirements for this AOAC validation study for the detection of penicillin G residues in bovine muscle tissue. The results also show that there is cross reactivity with many other antimicrobial compounds. PremiTest is a broad-spectrum screening test which can detect many relevant antimicrobials in meat.

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R-biopharm Premi®Test, AOAC Performance Tested MethodsSM Certification Number 060601

Table 1: Summary of performance parameters of Premi®Test for penicillin G obtained in the study (1)

Laboratory	Study	Performance Pen G (ppb)	Remarks
CSL	Part 1: 6 out of 6	10	A range of 5 different concentrations from 0 – 50 ppb was tested. At 10 ppb, all 6 of 6 samples analysed tested positive.
USDA	Part 2: 30 out of 30	10	Results confirmed that the detection sensitivity of the PremiTest for pen G was 10 ppb with 95% confidence level
DSM	Part 3: Ruggedness, interference, cross-reactivity	10	The presence of other antimicrobial drugs did not interfere with the PremiTest kits ability to detect pen G residues at a concentration of 10 ppb added to bovine muscle tissue.
USDA	Part 4: Incurred tissue	10	The detection sensitivity determined using fortified tissue samples was further confirmed using incurred bovine muscle tissues. PremiTest was able to detect all pen G residues in incurred samples containing pen G at concentrations greater than or equal to 12 ppb .

Table 9: Scanner results to confirm the Premi®Test's detection sensitivity as 10 ppb for penicillin G (1)

Zero control			
Sample		Sample	
	Interpreted Scanner Results		Interpreted Scanner Results
	NEG		
1	NEG	9	
2	NEG	10	
3	NEG	11	NEG
4	NEG	12	NEG
5	NEG	13	NEG
6	POS	14	NEG
7	NEG	15	NEG
8			NEG
Penicillin G 10 ppb			
Sample	Interpreted Scanner Results	Sample	Interpreted Scanner results
	POS		POS
	POS		POS
1	POS	16	POS
2	POS	17	POS
3	POS	18	POS
4	POS	19	POS
5	POS	20	POS
6	POS	21	POS
7	POS	22	POS
8	POS	23	POS
9	POS	24	POS
10	POS	25	POS
11	POS	26	POS
12	POS	27	POS
13	POS	28	POS
14		29	
15		30	

REFERENCE CITED

1. de Rijk, Angelique., Validation Study Report Premi®Test, AOAC Performance Tested MethodsSM certification number 060601.

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● Premi®Test detection limits in different animal food products

Substance	Chicken	Pork	Beef	Egg	Shrimp
β-lactams					
Amoxicillin	5	5	5	5	15
Ampicillin	5	5	5	5	
Penicillin-G	2,5	2,5	2,5	2,5	5
Cloxacillin		>100		100	
Oxacillin		100			
Dicloxacillin					
Cephalosporins					
Cefquinome	75	100	100		
Ceftiofur	100	200	100	400	
Macrolides					
Tylosin	50	25 - 50	50	50	
Erythromycin	100	100	100	50	100
Lincomycin	100	100	100		
Tilmicosin	50	50	50		
Spiramycin	1000	1000	1000		
Tetracyclines					
Chlortetracycline	100	100	100	600	1000
Oxytetracycline	100	100	100	400	100
Doxycycline	100	100	100	200	
Tetracycline		50		200	
Demeclocycline		50			
Sulphonamides					
Sulphamethazine	75	50-100	100	25	
Sulphadiazine	75	50-75	75	25	50
Sulphamethizole		50-100			
Sulphaguanidine	<200	150	<200		
Sulphadime-thoxine		25 - 50	<100		50
Sulphapyridine	<50	50	<100		
Sulphamethoxy-pyridine	<100	25			
Sulphisoxazole	<100	25			
Sulphathiazole	<100	25			
Sulphachloropyri-dazine	<100	25			
Sulphamerazine	<100	25	<100		
Sulphanilamide	<100	150			
Sulphaquinoxaline	<100	50	<50		
Sulphamethiazole	<100		<50		
Sulfamethoxazole				25	
Aminoglycosides					
Gentamicin	100	100	100	100	
Streptomycin	1500	1500	3000	1000	
Neomycine	300	300	300	300	200
Dihydrostrepto-mycin					
Quinolones					
Oxolinic acid					>10000
Enrofloxacin	>600	>600	>600		
Flumequine	>100	>100	>100		
Polypeptides					
Virginiamycin	500	500	500		
Bacitracin	500	500	500		
Zn-bacitracin	1250				
Colistin	>1000				
Ionophores					
Salinomycin	1000				
Monensin	1250				
Oligosaccharides					
Avilamycin	>5000				
Others					
Florphenicol	100	100	100		5000
Chloramphenicol	2500	2500	2500	2500	
Trimethoprim	50				
Narasin	1250				
Amprolium	>2000				
Phosphomycine	>1500				
Ronidazole					>5000
Furazolidone	>1500				

All detection limits are given in µg/kg = ppb

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Antibiotic residue screening

Premi®Test – microbial inhibition test for the screening of food

For food of animal origin

Simple to perform

Results in less than 4 hours

PERFORMANCE TESTED
Bovine muscle tissue was certified by the AOAC Research Institute.

RBP 31/02 - 4/11
ALTERNATIVE ANALYTICAL METHODS FOR AGRIBUSINESS
Certified by AFNOR Certification
End of validity: see AFNOR website
<http://nf-validation.afnor.org/en>
search for PremiTest

06/2026

Premi®Test – screening of antibiotic residues in animal origin food

In livestock breeding, antibiotics are used for the treatment of bacterial infections and diseases. Additionally, antibiotics can illegally be used as antimicrobial growth promoters (AGP). As a consequence, antibiotic residues can remain in food of animal origin and bear a health risk for the consumers.

The Premi®Test antibiotic residue screening test detects a broad spectrum of the most veterinary used antibiotics in animal husbandry. The easy-to-implement and easy-to-use test enables on-site screening by farmers, slaughterhouses, meat-processors, fisheries etc. In contrast to conventional plate tests, reliable results are available in less than four hours.

Premi®Test is suitable especially for the analysis of fresh meat.

Recommendations for sample-preparations of fish, shrimps, eggs, kidney, liver, urine, blood and feed are available. These matrices require individual validation by the customer.

Supportive information are available at R-Biopharm upon request.

Test kits and equipment

Premi®Test (Art. No. R3900): 100 test ampoules
or
Premi®Test (Art. No. R3925): 25 test ampoules

Premi®Test Starter Kit (Art. No. ZPT-2000)
Contains meatpress, incubator, timer, scissors
and briefcase

How to use Premi®Test – detection of antibiotic residues in meat

1 Remove the required amount of ampoules from the frame. Be careful not to damage the foil of the remaining ampoules. Use one additional ampoule for the negative control, which is a sample containing no antibiotics.

6 Incubate ampoules for 20 minutes at room temperature.

Switch on the Premi®Test Incubator for temperature stabilization at 64 °C.

2 Cut the meat into pieces of approx. 2 cm² and put one piece into the meatpress. Increase pressure slowly and hold constant until approx. 200 µL of meat juice is obtained.

7 Wash the ampoules twice with demineralized water. Remove excessive water by turning the ampoules upside down on a piece of soaking paper, but do not tap the ampoules.

3 Open the ampoules with one corner of the green frame and label the ampoules.

8 Close the ampoules with the supplied perforated foil.

4 Dip the pipette into the liquid and decrease the pressure on the bubble to allow the stem to fill with the meat juice.

Add the meat juice to the pipette by squeezing the smaller upper bulb.

9 Put the ampoules into the incubator at 64 °C.

5 By fully squeezing the top bulb, approximately 100 µL of the meat juice will be released.
(Note: A small excess of liquid may remain in the lower bulb.)

10 Read the results of the samples, when the negative control shows a clear colour change into yellow (approx. after 3 h).

Determination of results

Compare the colour of the sample ampoules to the colour of the negative control ampoule. If the sample ampoule remain purple or if the colour deviate from the colour of the negative control ampoule, the sample potentially contains antibiotic residues and the results have to be verified by the official reference method.

Premi®Test products

Premi®Test detection limits in different animal food products

Substances	Chicken	Pork	Beef	Substances	Chicken	Pork	Beef
β-lactams				Aminoglycosides			
Amoxicillin	5	5	5	Gentamicin	100	100	100
Ampicillin	5	5	5	Streptomycin	1500	1500	3000
Penicillin-G	2.5	2.5	10*	Neomycine	300	300	300
Cloxacillin		>100		Spectinomycine			5000
Oxacillin		100		Chinolones			
Dicloxacillin				Oxolin acid			
Cephalosporins				Enrofloxacin	>600	>600	>600
Cefquinome	75	100	100	Flumequine	>100	>100	>100
Ceftiofur	100	200	100	Polypeptides			
Macrolides				Virginiamycin	500	500	500
Tylosin	90	90	90	Bacitracin	500	500	500
Erythromycin	200	200	200	Zn-bacitracin	1250		
Lincomycin	>100	>100	>100	Colistin	>1000		
Tilmicosin	50	50	50	Ionophores			
Spiramycin	1000	1000	1000	Salinomycin	1000		
Tetracyclins				Monensin	1250		
Chlortetracycline	160	160	160	Lasalocid	10000		
Oxytetracycline	160	160	160	Oligosaccharides			
Doxycycline	100	100	100	Avilamycin	>5000		
Tetracycline		50		Others			
Demeclocycline		50		Florfenicol	100	100	100
Sulfonamides				Chloramphenicol	2500	2500	2500
Sulfamethazine	75	50-100	100	Trimethoprim	50		
Sulfadiazine	90	90	90	Narasin	1250		
Sulfamethizole		50-100		Amprolium	>2000		
Sulfguanidine	<200	150	<200	Phosphomycine	>1500		
Sulfadimethoxine		25 - 50	<100	Ronidazole			
Sulfapyridine	<50	50	<100	Furazolidone	>1500		
Sulfamethoxypyridine	<100	25					
Sulfisoxazole	<100	25					
Sulfathiazole	<100	25					
Sulfachloropyridazine	<100	25					
Sulfamerazine	<100	25	<100				
Sulfanilamide	<100	150					
Sulfaquinoxaline	<100	50	<50				
Sulfametiozole	<100		<50				
Sulfamethoxazole							

All detection limits are given in µg/kg = ppb.

* = Certified by the AOAC PTM Program.